# Rocket HSG & Sonohystogram Catheter INSTRUCTIONS FOR USE

**Indications:** For the instillation of all aqueous based x-ray and ultrasound contrast media to aid the diagnosis of Fallopian tube blockage, uterine abnormality and other pathological conditions associated with primary and secondary infertility. This device is only for use by or under the supervision of trained personnel and in accordance with local and national guidelines such as those published by the Royal College of Obstetricians & Gynaecologists and the Royal College of Radiologists.

**Contraindications:** Not for use in patients with a vaginal infection, sexually transmitted disease or suspected pregnancy.

**Preparation:**

- Using aseptic technique, open the package and remove the device and syringe. Inject 2.0ml of air into the balloon and lightly palpate with the fingers to confirm balloon patency. Withdraw **all** residual air from the balloon.
- This diagnostic procedure can normally be carried out without medication, although light IV sedation and/or analgesia may be necessary in very anxious patients.
- Place the patient into the lithotomy position with knees raised, expose the cervix with a Cusco or Sims type speculum. Assess the orientation of the cervix and establish the position of the external cervical os. Swab away any excess mucus with sterile normal saline or an aqueous antiseptic solution
- To aid insertion, a small quantity of water soluble lubricant may be applied to the tip and balloon.
- Pre-form the catheter with the wire stylet in place (applies only to R57422-00-ST and R57424-00-ST versions)
- Locate the distal tip of the Rocket® HSG catheter into the external cervical os. Using a very gentle rotating motion, advance the catheter into the uterus. The patient may report slight discomfort when the catheter passes through the internal os and into the uterine cavity. The outer sliding sheath may be used to provide extra rigidity to the catheter to aid initial insertion.
- Pre-fill the catheter with saline or contrast media to avoid the risk of erroneous diagnosis caused by air bubbles
- Attach the supplied pre-set syringe to the balloon inflation port and inflate with up to **2.0ml** of air or normal saline sufficient to create a seal (Fig. 1). Gently withdraw the catheter until it is located against the internal os. The pressure of the balloon against in the internal os creates a seal and prevents leakage of the contrast media.
- The catheter may also be inflated in the endocervical canal if insertion beyond the internal os is deemed undesirable.
- Attach a 10 or 20ml syringe to the rear port, ensuring the tap is in-line with the catheter. Instil 10-20ml of contrast media, sufficient to adequately demonstrate the uterine cavity and Fallopian tubes. Further volumes of contrast media may be added by turning the tap through 90° to the 'Off' position and replenishing the syringe.


Figure 1

**WARNING: Maximum balloon volume = 3.0 ml.**

Do not over inflate as this may cause severe discomfort or balloon rupture and loss of the contrast seal may occur.

Only suitable for use with aqueous based contrast medias

**CAUTION:** NOT for use with pressure injectors.

**CAUTION:** Balloon should be protected from light until use.

**To remove the device:** Attach the supplied syringe to the balloon port, open the tap and withdraw the plunger to deflate the balloon. Gently withdraw the catheter from the uterus.

**Disposal:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, 'those pertaining to human health & safety and care of the environment

**DO NOT RESTERILISE**

Unless opened or damaged, contents of package are sterile.

Store at room temperature. Avoid prolonged exposure to elevated temperatures.

This device is not manufactured with natural rubber latex

**Rx ONLY**

Manufactured in the UK by: **ROCKET MEDICAL PLC** Sedling Road, Washington, Tyne & Wear, NE38 9BZ England  
Rocket is a registered trademark of **ROCKET MEDICAL PLC.**

Distributed in the USA by **ROCKET MEDICAL**  
150 Recreation Park Drive Unit 3 Hingham MA 02043 USA Tel: 1-800-707-7625

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ZDOCK049 160517 Rev.15 Copyright© 1998-2017 **ROCKET MEDICAL PLC**