
# Rocket® Multipurpose Needle Drain

## INSTRUCTIONS FOR USE


**Scope:** These instructions cover all R58800-08-ND and R58800-08-NK Rocket Multipurpose Needle Drains, Kits and derivatives.

**Device Description:** The device comprises an 8FG x 20cm x 6 hole polyurethane pigtail catheter with flexible support and stitch plate loaded over a 17G cannula/needle assembly.

**Indications:** For external drainage for ascites or abdominal cavity abscess by one-step puncture. This device should only be used by, or under the supervision of, appropriately trained personnel in conjunction with national clinical practice guidelines such as those published by the Royal College of Radiologists and Royal College of Surgeons.

**Contraindications:** Not to be used where the risks of insertion outweigh the benefits of drainage or where the anatomy cannot be adequately determined by ultrasound or fluoroscopy.

### Procedure:

- Following local hospital guidelines, use aseptic technique to prepare the site for insertion.
- For drainage of abdominal ascites, use the diagram (right) as a guide to insertion sites.
- Prepare the drain for insertion by sliding the red sleeve forward to straighten the pigtail.**
- Using a gentle rotating action, advance the needle set forward through the pigtail.
- When fully inserted, lightly push the hubs together, stretching the catheter tube until the needle bevel tip is seen protruding through the catheter tip.
- Lock the needle set to the flexible support with a clockwise twisting action.
- REMOVE THE RED STRAIGHTENING SLEEVE and discard before insertion.**
- If the drain is to be inserted under aspiration, remove the rear cap and attach a 10-20ml syringe to the needle hub.
- Having chosen the site, it is strongly recommended that a 2-3mm skin nick is used prior to insertion.
- Grasp the scalpel and carefully extend the blade by moving the slide toward the tip of the scalpel, using the thumb of the hand holding the scalpel. Extend the slider until you reach the positive stop, the slider will locate into a notch when it is completely extended. To retract the blade, grasp the scalpel carefully and move the slide toward the back of the scalpel, using the hand holding the scalpel. You should feel clicks as the blade is retracted and a positive stop once the blade is completely retracted.
- To retract the blade permanently, move the slider past the notch at the back of the scalpel.
- It is recommend that the needle drain is gripped by the stitch plate and flexible stem only.
- Advance the drain only until a loss of resistance is felt or aspiration confirms insertion into the cavity.
- Unlock the needle set with an anti-clockwise twisting action and withdraw the needle set *whilst advancing the catheter*.
- Do NOT remove the needle set until correct insertion is established. Once correct insertion is established the drain may be fully advanced into the cavity
- Do NOT attempt to re-insert the needle set – remove the whole drain and start again from #3.**
- The position of the needle point relative to the catheter end is detailed in the table above. The needle point is *inside* the catheter tip when the locking ring clears the luer connector
- Connect the drain connector set and attach to an appropriate collection bottle or bag.
- Confirm the catheter is draining freely before completing the procedure.
- The device can be secured either by suturing the stitch plate into position or by applying dressing strips over the flats of the stitch plate and securely fixing to the skin.

Needle point position	Hub locking ring position

**WARNING: Do NOT attempt to insert the needle into the drain without first straightening the pigtail with the protective sleeve.**

**If the needle is allowed to puncture the catheter wall, the drain should be discarded and replaced with a new one.**

**CAUTION: Where anatomy is potentially distorted or compromised, insertion should be conducted under continuous imaging such as ultrasound or radiological screening.**

**This product has been coated with silicone and beads may form on the needle surface during storage.**

**Disposal:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

Rx ONLY

### NOT FOR LONG TERM USE.

Unless opened or damaged, contents of package are sterile.

Store at room temperature. Avoid prolonged exposure to elevated temperatures.

Manufactured in the UK by: **ROCKET MEDICAL PLC**. Sedling Road. Washington. Tyne & Wear. England. [www.rocketmedical.com](http://www.rocketmedical.com)

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Copyright© 2007-17. **ROCKET MEDICAL PLC**
All Rights Reserved. (GB)