# Rocket® FBS Kits with Sampling Wand

## INSTRUCTIONS FOR USE

**Indications: Fetal Blood Sampling (FBS) Kit.** Procedure kit for single use, to obtain up to fetal blood pH/gas samples in cases of suspected fetal distress during labour. If additional sampling is necessary or if labour is protracted, a further pack may be required. This device should only be used by or under the supervision of trained personnel and in conjunction with current clinical practice guidelines such as: **NICE Guidelines No CG190:2017 'Intrapartum care of healthy women and their babies during childbirth' (as amended).**


These instructions apply to all Rocket® Fetal Blood Sampling Kits with Sampling Wands. Contents of packs will vary. Rocket Medical supply a range of Rocket FBS Kits suitable for a wide range of blood gas analysers, contact Customer Services for advice on compatibility.

**Contraindications:** This procedure is contraindicated with placenta previa, when identification of the presenting part is uncertain, in the presence of genital infections (e.g. herpes, Group B streptococcus, gonorrhoea, HIV/AIDS and Hepatitis B) or where the mother is a confirmed carrier of haemophilia and the fetus is either affected or of unknown status.


FBS Kit	Description	Volume $\mu$ l	Tube Length	OD	ID	Heparin type	Concentration usp/ml
R57012	Straight Glass, Heparinised	2320	295 $\pm$ 5.0mm	3.0 $\pm$ 0.15mm	1.0 $\pm$ 0.04mm	Na	100
R57013	Angled Glass, Heparinised	2320	295 $\pm$ 5.0mm	3.0 $\pm$ 0.15mm	1.0 $\pm$ 0.04mm	Na	100
R57012/013	Plastic Capillary, Heparinised	75	74.85mm	1.87mm	1.08mm	Ca/Na/Li	75
R57018/028	Plastic Capillary, Heparinised	55	59mm +2.0 - 0.0mm	2.05 $\pm$ 0.1mm	1.25 $\pm$ 0.05mm	Na/Li	240
R57024	Plastic Capillary, Heparinised	85	87.0 $\pm$ 0.15mm	2.02 $\pm$ 0.01mm	1.17 $\pm$ 0.01mm	Na/Li	240
R57026	Plastic Capillary, Heparinised	100	85.0mm	2.3 $\pm$ 0.05mm	1.20 $\pm$ 0.05	Ca/Na/Li	200
R57027	Plastic Capillary, Heparinised	170	100.0mm	2.0 $\pm$ 0.05mm	1.55 $\pm$ 0.05	Ca/Na/Li	100

### Procedure:

- Using aseptic technique, remove the kit from outer bag, place onto a suitable surface and unfold outer drape to form a sterile field.
- Select a capillary tube from the tube in the kit and load into the holder. Locate the end of the tube into the holder at ❶ and then gently advance into the rear tube grip at ❷. Ensure the tube is secure before sampling.
- Follow local hospital procedure, prep the perineum and vagina and complete a vaginal examination to establish the cervical dilatation and station of the fetal head.
- Insert the plastic capillary tube inserted into the sampling wand tube holder, apply a *light* pressure to secure the capillary tube. The kit contains a secondary double ended capillary tube holder if preferred.
- Locate the Amnilume™ light module into the slot in the amnioscope and push gently to engage. The light will activate when the light source is correctly fitted. To turn the light off, withdraw the light module.
- Insert amnioscope complete with obturator into the vagina against the fetal presentation ensuring that the cervix is not trapped.
- Remove obturator. Position the bevelled end by rotating endoscope to obtain the best possible seal against the fetal presentation. Rotate to finally position the light source module in the upper quadrant.
- Maintain close contact between amnioscope and scalp to prevent leakage of liquor.
- Clean puncture area with the green x-ray detectable gauze swabs and forceps provided.
- Spray the fetal scalp with ethyl chloride spray or other suitable agent to promote capillary action if required.
- With a clean swab stick, apply a *thin* layer of white soft paraffin (supplied) to puncture area to promote capillary blood aggregation. **Remove excess with a new dry swab stick before skin puncture.**

- The sampling wand comes pre-fitted with an FBS blade. Lightly holding the shaft and rear finger grip (A), push the grip forward to expose the blade (B). **Ensure the blade is locked forward before use.**
- Perform skin puncture through the amnioscope. Either ❶ make a vertical stab, rotate the blade 90° and make a second stab to form a cross or ❷ make a vertical stab, then angle the blade downwards and extend the incision using the front edge of the blade to produce a single incision 4-5mm long.
- Always wait for a good sized (3-4mm) bead to form on the surface before collection. Never slash at the scalp and avoid multiple small stabs.
- To collect the blood sample, depress the button on the wand shaft to retract the FBS blade and angle the capillary tube into the bead.
- Rapidly collect blood sample, filling the tube completely, **without air bubbles.**
- Twist lightly to remove the capillary tube from sampling wand and cap both ends with covers supplied. Take to the analyser immediately.
- Do NOT attempt to mix, shake, tip or otherwise disturb the sample. Rocket® FBS Kits contain high heparin tubes, mixing is not required or recommended
- Always ensure the **clean**, non-fetal exposed end is presented to the analyser. Please refer to the analyser instructions for guidance.
- Always ensure that the **correct** measuring mode is selected *before* performing analysis.

**Post sampling:**

21. Haemostasis is accomplished by applying pressure to puncture site with a swab for one full contraction.
22. If test results so indicate, repeat procedure. The kit is designed to take up to 5 individual samples on the same patient.

**CAUTION: Avoid fetal face, fontanels or genitalia.** Rocket® FBS blades have a maximum depth of penetration of 2.2mm

**CAUTION: Do not attempt to shake, roll or tilt the filled collection tube or use 'fleas' and magnets to mix the heparin. Use of these techniques may cause air entrainment into the sample. Samples containing air bubbles may be rejected by the analyser.**

**CAUTION: Failure to use correct blood volume, select correct analysis mode or the use of clotted samples or those contaminated with air bubbles may lead to errors in results or rejection of the whole sample. Carefully follow the analyser manufacturers instructions.**

**Disposal:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment. The Rocket® Amnilume™ light module contains lithium-ion based power cells.

This device is not manufactured with natural rubber latex

**R<sub>x</sub>**  
ONLY

Li-ion

---

Unless opened or damaged, contents of package are sterile.  
Store at room temperature. Avoid prolonged exposure to elevated temperatures.

---

Manufactured in the UK by: **ROCKET MEDICAL PLC** Sedling Road, Washington, England, NE38 9BZ  
Rocket® is a registered trademark of **ROCKET MEDICAL PLC**.

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Copyright© 2010-17. **ROCKET MEDICAL PLC**