Rocket® Fetal Blood Sampling Kits

INSTRUCTIONS FOR USE

Indications: These instructions apply to Rocket Basic Fetal Blood Sampling Kits and derivatives. Contents of packs will vary. Procedure kit for single use, to obtain up to 5 fetal blood pH/gas samples in cases of suspected fetal distress during labour. If additional sampling is necessary or if labour is protracted, a further pack may be required. This device should only be used by or under the supervision of trained personnel and in conjunction with current clinical practice guidelines such as: **NICE Guidelines No. CG190: 02/2017.**

'Intrapartum care of healthy women and their babies during childbirth' (as amended).

Contraindications: This procedure is contraindicated with placenta previa, when identification of the presenting part is uncertain, in the presence of genital infections (e.g. herpes, Group B streptococcus, gonorrhoea, HIV/AIDS and Hepatitis B), where the mother is a confirmed carrier of haemophilia and the fetus is either affected or of unknown status or acute fetal compromise with associated fetal bradycardia of >3 mins.

FBS Kit	Description	Volume μ l	Tube Length	OD	ID	Heparin type	Concentration usp/ml
R57012	Straight Glass, Heparinised	2320	295 \pm 5.0mm	3.0 \pm 0.15mm	1.0 \pm 0.04mm	Na	100
R57013	Angled Glass, Heparinised	2320	295 \pm 5.0mm	3.0 \pm 0.15mm	1.0 \pm 0.04mm	Na	100
R57012/013	Plastic Capillary, Heparinised	75	74.85mm	1.87mm	1.08mm	Ca/Na/Li	75
R57018	Plastic Capillary, Heparinised	55	59mm +2.0 - 0.0mm	2.05 \pm 0.1mm	1.25 \pm 0.05mm	Na/Li	240
R57024	Plastic Capillary, Heparinised	85	87.0 \pm 0.15mm	2.02 \pm 0.01mm	1.17 \pm 0.01mm	Na/Li	240
R57026	Plastic Capillary, Heparinised	100	85.0mm	2.3 \pm 0.05mm	1.20 \pm 0.05	Ca/Na/Li	200
R57027	Plastic Capillary, Heparinised	170	100.0mm	2.0 \pm 0.05mm	1.55 \pm 0.05	Ca/Na/Li	100

Rocket Medical supply a range of Rocket® FBS Kits suitable for a range of blood gas analysers; contact Customer Services for advice on compatibility.

Preparation:

- Using aseptic technique, remove the kit from outer bag, place onto a suitable surface and unfold outer drape to form a sterile field.
- Follow local hospital procedure, prep the perineum and vagina and complete a vaginal examination to establish the cervical dilatation and station of the fetal head.
- Unpack one FBS blade and using a clockwise rotating action, locate the screw securely into the blade holder.
- Locate the Amnillum™ light module into the slot in the amnioscope and push gently to engage. The light will activate when the light source is correctly fitted. To turn the light off, withdraw the light module.
- Insert amnioscope complete with obturator into the vagina against the fetal presentation ensuring that the cervix is not trapped.
- Remove obturator. Position the bevelled end by rotating endoscope to obtain the best possible seal against the fetal presentation. Maintain close contact between amnioscope and scalp to prevent leakage of liquor
- Clean puncture area with the green x-ray detectable gauze swabs and forceps provided.
- With a clean swab stick, apply a *thin* layer of petroleum jelly (supplied) to the puncture area to promote blood aggregation. **Always remove excess petroleum jelly with one or two new, dry swab sticks before performing the skin puncture.**
- Perform skin puncture through the endoscope using the loaded FBS blade holder. Either ① make a vertical stab, rotate the blade 90° and make a second stab to form a cross or ② make a vertical stab, then angle the blade upwards and extend using the front edge of the blade to produce a single incision 4-5mm long. Always wait for a good sized (3-4mm) bead to form on the surface before collection. Never slash at the scalp and avoid multiple small stabs.

CAUTION: DO NOT USE ETHYL CHLORIDE SPRAY – this may adversely affect the internal surface of the plastic amnioscope

CAUTION: Avoid fetal face, fontanels or genitalia. Rocket® FBS blades have a maximum depth of penetration of 2.2mm

- Where supplied and approved by the analyser manufacturers, use a long glass heparinised tube to collect the bead of aggregated blood. Collection needs to be expedited rapidly to avoid blood clotting. **CAUTION: Do NOT use the sample if the blood clots**
- For micro volume analysers requiring precision volume samples of 170 μ l or less, use the correct micro-volume capillary tube inserted into the capillary tube holder. Rapidly collect blood sample, filling the tube completely, **without air bubbles.**
- Lightly twist the capillary to remove it from holder and cap both ends with covers supplied. Take to the analyser immediately.

CAUTION: Do not shake the filled collection tube or use wires and magnets to mix heparin.
Use of these techniques may cause air entrainment into the sample. Samples containing air bubbles may be rejected by the analyser.

- Always ensure the **clean**, non-fetal exposed end is presented to the analyser. Please refer to the analyser instructions for guidance. Always ensure that the **correct** measuring mode is selected *before* performing analysis.

CAUTION: Failure to use correct blood volume, select correct analysis mode or the use of clotted or contaminated samples may lead to errors in results or rejection of the whole sample. Carefully follow the analyser manufacturer's instructions.

Post sampling:

- Haemostasis is accomplished by applying pressure to fetal puncture site with a swab for one full contraction.
- If test results so indicate, repeat procedure.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment. The Rocket® Amnilume™ light module contains lithium-ion based power cells.

This device is not manufactured with natural rubber latex

Rx ONLY

STERILE EO

Li-ion

ROCKET MEDICAL PLC Sedling Road, Washington,
England, NE38 9BZ www.rocketmedical.com

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

DO NOT RESTERILISE

Unless opened or damaged, contents of package are sterile.
Store at room temperature. Avoid prolonged exposure to elevated temperatures.

0088

ZDOCK066

Rev. 30 300719

Copyright© 2010-19.

ROCKET MEDICAL PLC

All rights reserved (GB)

“Uncontrolled Copy – subject to change at any time without notification”

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Copyright© 2010-19. **ROCKET MEDICAL PLC**
Copyright© 2010-19. ROCKET MEDICAL PLC Wszelkie prawa zastrzeżone.