
Product Name: Water Trap Set for CRAFT™ Suction Pumps

ZDOCK150 Rev. 12

Language	Page
EN	2-4
FR	5-7
NL	8-10
IT	11-13
CN	14-16

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Australian Sponsor:

ROCKET MEDICAL PTY LTD
Suite 209, 1 Katherine Street,
Chatswood, NSW 2067 Australia

Rocket, ROCKET, the R device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc.

Water Trap Set for CRAFT™ Suction Pumps INSTRUCTIONS FOR USE

Product Name: Water Trap Set for CRAFT™ Suction Pumps

Product Code: R57685

Device Image:

Intended Use: Intended for the connection of oocyte aspiration needles to the CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps.

Indications: For the connection of oocyte aspiration needles to the CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps permitting vacuums up to 400 mmHg for the recovery of oocytes as part of the treatment of infertility. This device may only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines.

Contraindications: Not intended to be used in conjunction with any other oocyte aspiration pumps except Rocket Medical CRAFT™ and CRAFT™ DUO-VAC™ Pumps.

WARNING: Do not allow set to become contaminated with fluid. If set becomes contaminated replace immediately.

WARNING: Vacuum settings should always be made with the patient connection tube fully occluded.

Precaution: **Potential for DEHP Exposure from use of this device.** DEHP may leach from this device, which may cause reproductive harm in male neonates, pregnant women carrying male fetuses, and peripubertal males. The extent of exposure largely depends upon the medical treatments administered and the duration of the treatment. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure. (For further details please see 'POL58' at www.rocketmedical.com.)

Benefit Risk: Rocket Medical Plc has taken all necessary steps to ensure that residual risks associated with the use of Rocket Water Trap Set are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical Plc concludes that the overall medical benefits of the Rocket Water Trap Set, outweigh the possible risks when used according to the intended use.

Performance Claims and Clinical Benefit:

Performance Claim: To protect CRAFT™ and CRAFT™ DUO-VAC™ Suction pumps from the ingress of fluid

Clinical Benefit: Indirect clinical benefit, when used in conjunction with aspiration needles and CRAFT™ or CRAFT™ DUO-VAC™ Suction pumps for the harvesting of oocytes for Assisted Reproductive Technology procedures

Undesirable Side Effects: Bleeding, pain, infection and haematuria

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Description of the device: This device comprises of a clear glass bottle (to enable detection of aspirate should it enter the trap) and a tubing set. The patient tube set is 2m in length and has an integral Luer lock connector to ensure secure, reliable connection to an aspiration needle set.

User: A healthcare professional proficient in the undertaking of oocyte aspiration, working in accordance with local and national guidelines. The device may also be used by a trainee under the supervision of a healthcare professional proficient in the undertaking oocyte aspiration.

Patient Population: Adult patients undergoing oocyte aspiration for the treatment of infertility with IVF or other gynaecological patients such as those undergoing oocyte retrieval for cryopreservation or those with ovarian hyperstimulation syndrome.


Environment: To be uses in a sterile clinical setting.

For Single Use Only: the device may not have structural integrity to support reuse, resulting in harm to the patient and/or user. Re-use may result in harm to patients and users from cross-contamination and infection.

Sterile: the device is supplied sterile by Ethylene Oxide (EO) . Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: This device is not intended to be routinely found in an MRI environment and compatibility has not been verified.

Symbols and safety signs used on the device and labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Do not re-use
	Date of manufacture		Consult the instructions for use
	Use-by date		Caution
	Batch code		Does not contain natural rubber latex
	Catalogue number		Contains or presence of phthalate
	Sterilized using ethylene oxide		Medical device
	Do not resterilise		Do not use if package is damaged and consult instructions for use

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

	Single sterile barrier system with protective packaging outside		Single sterile barrier system
	Distributor		Importer

Operating Instructions:

CHECK: do not use if the packaging is open or damaged or the device is damaged.

PREPARE: Remove and discard reusable bottle set if present, including bottle, grey stopper, metal tubes, silicone tubing and adapter (applies to older models only).

USE:

1. Remove and discard the reusable bottle set if present, including bottle, grey stopper, metal tubes, silicone tubing and adapter (applies to older models only).
2. Observing aseptic techniques, open the Water Trap Set for CRAFT™ Suction Pumps onto a sterile field and pass the bottle out to the circulating assistant. Retain the distal end of the long patient connection tube in the sterile field.
3. Carefully insert the bottle into the bottle holder (see figure 2).
4. Locate the short length of tube and attach firmly to the suction inlet (indicated in figure 3). Ensure there are no kinks or restrictions in the tube.
5. Fully occlude the *patient* connection tube between forefinger and thumb downstream of the water trap bottle.
6. Ensure the pump is connected to the correct power source as indicated in ZDOCK105.
7. Maintain the occlusion and operate the Standard Vacuum footswitch.
8. Adjust the vacuum control to set the vacuum to between 100-200mmHg.
9. *Note: Collection needles with tube sets >70cm in length typically require vacuums >150mmHg to achieve satisfactory flow rates.*
10. Once set, on CRAFT™ DUO-VAC™ models, operate the high vacuum footswitch, or on CRAFT™ Suction models, occlude the red button (see figure 4), to achieve maximum vacuum. The vacuum reading should exceed 400mmHg⁻¹.
11. Release the tube set from between finger and thumb, The system is now ready for use.
12. Do NOT allow the bottle to fill with fluid. Replace immediately if the bottle becomes contaminated.

Figure 2: Bottle Holder

Figure 3: Suction Inlet

Figure 4: Red Button on CRAFT™ DUO-VAC™

Disposal: This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Ensure sharps are disposed of in sharps containers.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pnf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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