Product Name: Copeland® Fetal Scalp Electrode
ZDOCK167 Rev. 15

Language	Page
EN	2-5
DE	6-9
NL	10-13
SV	14-17
FR	18-21
ES	22-25

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Australian Sponsor:

ROCKET MEDICAL PTY LTD
Suite 209, 1 Katherine Street,
Chatswood, NSW 2067. Australia

Copeland, the R device and the Rocket Medical logo are registered trademarks of Rocket Medical Plc. Philips and Avalon are registered trademarks of Koninklijke Philips N.V.

Product Names:

Copeland® Fetal Scalp Electrode for Qwik Connect Cables

Copeland® Fetal Scalp Electrode for Philips® Avalon® FM30/FM50 DECG Cables

Product Codes:

R57008-00-CN

R57008-00-FM

Device Image:

Intended Use: Intended for attachment to the fetal scalp to enable transvaginal fetal heart rate monitoring.

Indications: Indicated for use during labour, where transvaginal fetal heart rate monitoring is required.

Contraindications:

Not indicated in the following situations:

- where the amniotic membranes are intact and artificial rupture is not indicated.
- in cases of major placenta previa as this may cause significant maternal bleeding.
- in cases of face presentation as there may be permanent facial scarring or injury to the eye.
- when identification of the presenting part is uncertain;
 - attachment to the scalp overlying the fontanelles may lead to loss of CSF if the meninges are penetrated.
 - attachment to the face may lead to damage to the fetal eye and cause blindness
 - attachment to the face may cause permanent scarring.
 - attachment to the fetal genitalia in an undiagnosed breech presentation may cause damage to the scrotum or labia.
- where there is a risk of vertical or perinatal transmission (for example, active herpes simplex virus, HIV and hepatitis C). Use caution and refer to your local clinical policies and procedures in cases when the woman is positive for Group B streptococcus, gonorrhoea, hepatitis A, hepatitis B or syphilis.
- where the mother is a confirmed carrier of haemophilia and the fetus is either affected or of unknown status as this may lead to increased fetal bleeding at the application site.

Not indicated when the fetus is premature (less than 34+0 weeks), unless all of the following apply:

- it is not possible to monitor the fetal heart rate using either external cardiotocography or intermittent auscultation
- it has been discussed with a senior obstetrician
- the benefits are likely to outweigh the potential risks
- the alternatives (immediate birth, intermittent ultrasound and no monitoring) have been discussed with the patient and are unacceptable.

Use in premature fetuses may lead to or worsen existing fetal anaemia if there is bleeding after application of the scalp electrode. In fetuses between 34+0 and 36+0 weeks of pregnancy, consideration should be given to external CTG or intermittent auscultation prior to a decision to use internal fetal monitoring.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

WARNING: Do not apply the electrode to the fetal face, fontanelles or the genitalia.

In the event of a poor trace or no trace at all, follow this sequence until a trace is established.

1. Check that you have followed the Instructions for Use and applied the fetal scalp electrode correctly.
2. Check that the fetal monitor and paper recorder is switched on.
3. Check that the maternal attachment electrode is correctly applied and is connected to the leg plate cable.
4. Check that the leg plate cable is correctly connected to the monitor. If using wireless transducers, check that there are no issues with the wireless connection.
5. Check that the fetal scalp electrode is firmly and securely inserted into the leg plate cable.
6. For Philips DECG cables, confirm that the cable port has been maintained in accordance with the manufacturer's cleaning instructions.
7. Re-position the electrode.
8. Replace the electrode.

WARNING: Chronic fetal bleeding (e.g., due to partial placental abruption) leading to loss of fetal blood volume, may compromise the fetus's ability to successfully respond to hypoxia.

CAUTION: Do not over-rotate or apply excessive torque when applying or removing the device. This may cause damage to the spring/hook assembly, resulting in poor scalp attachment or difficulty detaching the device.

CAUTION: Do not rely on the spring tension to insert the needle into the skin.

CAUTION: Device performance and safety has not been verified for use in water births. The device must not be immersed in water.

CAUTION: If there is severe maternal thrombocytopenia, any decision to use the fetal scalp electrode should be taken carefully, considering the potential risks and benefits, and taking account of clinical guidelines. In case of accidental maternal injury from the electrode, or if the fetus is affected, there may be increased bleeding, bruising or abnormal blood clotting.

Benefit Risk: The clinical benefit of the device (determination of the fetal heart rate to aid in decision making during labour) outweighs the possible risks when used according to the intended use.

Performance Claims and Clinical Benefit:

The intended performance claims for the device are:

- Attaches to the fetal scalp
- Remains attached to the fetal scalp
- Enables fetal heart rate monitoring

The intended clinical benefit is:

- Determination of the fetal heart rate to aid in decision making during labour.

Undesirable Fetal Side Effects: Bleeding, scalp trauma, scarring, hematoma, inflammation and infection.

Description of the Device: These devices enable the transmission of the fetal heart rate from the fetal scalp to a CTG system, for the monitoring of the fetal heart rate, as a component of CTG monitoring.

User: The intended user is a healthcare professional proficient in the undertaking of fetal heart rate monitoring or a person in training under supervision of a suitable competent healthcare professional, working in accordance with local and national guidelines.

Patient Population: The intended patient demographic is pregnant patients during labour, where internal monitoring of the fetal heart rate is required.

Environment: This device is intended to be used in a clinical setting. No other environmental requirements or restrictions apply.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

For Single Use Only: The device may not have structural integrity to support reuse (other than repositioning), resulting in harm to the patient and/or user. Re-use (other than repositioning) may result in harm to patients and users from cross-contamination and infection.

Sterile: The devices are supplied sterile by gamma irradiation. Do not re-sterilise; it may not be possible to achieve the required sterility assurance level upon re-sterilisation of a used device. Re-sterilisation may also compromise the structural integrity of the device, leading to device failure.

MRI Safety: MRI compatibility has not been verified.

General Information:

R57008-00-CN is compatible with Qwik Connect Cables.

R57008-00-FM is compatible with Philips® Avalon® FM30/FM50 DECG Cables.

The device has been assessed for use for periods not exceeding 24 hours.

Symbols and Safety Signs Used on the Device and Labelling:

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Indicates the Authorized Representative in the European Community		Keep dry
	Date of manufacture		Do not re-use
	Use-by date		Caution: Federal law restricts this device to sale by or on the order of a physician.
	Importer		Consult the instructions for use
	Batch code		Caution
	Catalogue number		Does not contain natural rubber latex
	Sterilized using irradiation		Medical device
	Do not re-sterilise		Do not use if package is damaged and consult instructions for use
	Distributor		Importer
	Single sterile barrier system		Single sterile barrier system with protective packaging outside

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Operating Instructions:

CHECK: Do not use if the packaging is open or damaged, or the device is damaged.

PREPARE:

1. Conduct a vaginal examination to confirm that it is appropriate to use the fetal scalp electrode and to identify the desired attachment site.
2. Advance the fetal scalp electrode between the examining fingers until the needle (A) is held against the fetus using light finger pressure on the serrated depression of the finger grip (B).

APPLICATION AND CONNECTION:

3. With the fingers of the other hand, rotate the finger control (C) approximately 180° ANTI-CLOCKWISE to open the needle and then approximately 180° CLOCKWISE to insert the needle into the fetal scalp. Take care not to use excessive force or to over-rotate the needle, as this may damage the mechanism, resulting in poor scalp attachment or difficulty detaching the device. Do not rely on the spring tension to insert the needle into the skin.
4. **Firmly** connect the fetal scalp electrode to the leg plate cable and the leg plate cable to the maternal attachment electrode.
5. Follow the monitor manufacturer's instructions for use/user manual to connect the leg plate cable to the fetal monitoring system.




REMOVAL:

6. To remove the fetal scalp electrode, turn the finger control (C) 180° ANTI-CLOCKWISE. It is not necessary to hold the finger grip (B).

REPOSITIONING:

7. The device may be repositioned, if required, by following the preparation and application steps above.

Disposal: This device should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment. Failure to do so may increase the risks of infection or other microbial hazards. Take care when handling sharps, to avoid needlestick injuries. Ensure sharps are disposed of in sharps containers.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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