INSTRUCTIONS FOR USE

ZDOCKI80

Rev. 21

	Page
EN	2-5
DE	6-9
FI	10-13
FR	14-17
GK	18-21
NL	22-25
NO	26-29
PT	30-33
SE	34-37
AR	38-41

Rocket Seldinger Chest Drain Kit

INSTRUCTIONS FOR USE

(EN)

Scope: These instructions cover the **R54544-12-SG, R54544-18-SG, R54544-12-PK, R54544-18-PK, R54549-16-PK, R54549-20-PK, R54549-24-PK, R51549-18-PK and R51549-20-PK** Rocket Seldinger Chest Drainage Sets.

Product Description: The Rocket Seldinger Chest Drainage Sets **R54544-12-SG and R54544-18-SG** contains a 10ml syringe, introducer needle, wire guide, No 11 safety Blade Scalpel, 10cm SG dilator(s) with initial 5cm active length, Luer fir tree and Luer wide bore tubing set connectors, Luer 3 way tap, coaxial chest tube inserter and clear drainage catheter with radiopaque stripe, needle safety device and catheter fixation device. In addition to the above the **R54544-12-PK, R54544-18-PK, R54549-16-PK, R54549-20-PK, R54549-24-PK, R51549-18-PK and R51549-20-PK** make up a complete procedure tray with a disposable gown, fenestrated drapes, skin prep sponges and gallipot, 20ml syringe, safety green needle and safety orange needle.

Indications: This product has been designed for the percutaneous introduction of a chest tube for drainage of air and or pleural fluid.

Contraindications: Non-pleural drainage.


The **12Fg & 16Fg** catheters are **not** suitable to be used for drainage of viscous fluids.

WARNING:

- **Over insertion of the dilator into the chest cavity is associated with serious injury and MUST be avoided.**
- **Only insert the dilator sufficient to dilate the superficial track.**
- **This device should only be used by, or under the supervision of, appropriately trained personnel and in conjunction with current local clinical practice guidelines**

Instructions for Use:


1. Following local hospital policy, and wearing the sterile gown, prepare the chest tube insertion site with an approved solution such as povidone-iodine and drape with fenestrated drape to maintain aseptic technique.
2. Identify and mark the insertion point following current practice guidelines. After injection of local anaesthetic, make a small skin incision 4-5mm.

CAUTION: Insertion should be just ABOVE the rib to avoid damage to the intercostal neurovascular bundle.


When removing the protective sheath on the access needle, hold the silicone ring to prevent it from being removed. The measurement provided by using the silicone ring can be used to determine if the correct length of dilator is being used, and if the safety guard needs to be removed.

3. Attach the introducer needle to a syringe filled with a local anaesthetic e.g. lignocaine, and administer additional local anaesthetic to the skin and underlying tissues.
4. Advance the introducer needle over the superior border of the rib and into the pleural space. Fluid or air should be aspirated to verify correct position.
5. The needle should be introduced and directed with appropriate orientation inferiorly or superiorly.
6. Remove the syringe and advance the soft “J” end of the wire guide through the needle and into the pleural space. Before removing the needle, move the silicone measurement ring down to the patient’s skin.

CAUTION: The wire guide should pass through the needle and advance into the pleural space without resistance.

7. Remove the introducer needle.
8. Using the measurement obtained using the depth marker on the introducer needle; measure against the safety guard dilator to ensure the correct length of dilator is being used. If the depth required is longer than the preset 50mm then remove the Safety Guard as described below.


Removing the safety guard:

CAUTION: Ensure sufficient length of wire guide remains external to the chest to facilitate controlled introduction of the dilators and drain.


- 9. While maintaining the wire guide position, dilate the track and opening into the pleural space by advancing and withdrawing, one by one, the supplied dilator(s) over the wire guide. Introduction into the pleural space is facilitated by rotating and advancing the dilators in the same plane as the wire guide to prevent kinking.
- 10. Remove the final dilator from the wire guide
- 11. Locate the drain set. Insert the free end of the guide wire into the distal catheter tip ensuring it travels through the stiffening tube.
- 12. With the wire guide still positioned within the pleural space, advance the chest tube inserter/chest tube assembly over the wire guide and into the pleural space.

CAUTION: It is important to advance the chest tube assembly into the pleural space whilst maintaining its direction in the same plane as the wire guide. This will make introduction easier and avoid kinking of the wire guide.

CAUTION: Ensure all side holes of the chest tube are positioned within the pleural space.

- 13. Remove the wire guide and chest tube inserter leaving the chest tube in place.
- 14. The chest tube can now be secured to the skin by using the fixation device R54960-PG-FX or suture or both and is ready for connection to a suitable sealed drainage system such as the Rocket Blue Bottle or Rocket Ambulatory Bag
- 15. Confirm the correct positioning of the drain with imaging as soon as possible.

How to use the Swann Morton Safety Scalpel

- 1. Grasp the scalpel and carefully extend the blade by moving the slide toward the tip of the scalpel, using the thumb of the hand holding the scalpel.
- 2. Extend the slider until you reach the positive stop, the slider will fit in the notch when it is completely extended.
- 3. To retract the blade, grasp the scalpel carefully and move the slide toward the back of the scalpel, using the hand holding the scalpel.
- 4. You should feel clicks as the blade is retracted and a positive stop once the blade is completely retracted.
- 5. To retract the blade permanently, move the slider past the notch at the back of the scalpel.
- 6. Dispose of the Swann Morton Retractable Disposable Scalpel in a puncture resistant container approved for sharps disposal, or in a manner consistent with your hospital procedure.

How to apply the Catheter Fixation Device

Indications: For use with the Rocket Seldinger catheter.

Contraindications: Not for use with catheters other than included in the Rocket Seldinger pack.


Device Description: The Fixation device can be used as secondary dressing, use of a suture to secure the catheter is recommended. This dressing is not designed to replace a suture to hold the catheter in place. The dressing will normally stay in-situ for 4 – 5 days, or up to seven days if required.

Procedure: For application of the drain fixation device follow the diagrams below.


The catheter fixation device enables the following:

1. Reduction of risk from an accidental displacement of catheter,
2. Reduction of risk of the catheter being expelled by the body
3. Safe anchoring and protection of catheter for up to one week.
4. The patient to shower with catheter inserted.
5. Hydrocolloid carrier absorbs exudate from puncture site
6. Prevents contamination of puncture site from environment.
7. Prevents kinking of catheter

Catheter fixation instructions for catheters 16Fg to 30Fg

Catheter fixation instructions for catheters 5Fg to 16Fg


CONTINUOUS USE SHOULD NOT EXCEED 28 DAYS

Disposal: This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

ROCKET MEDICAL PLC
Sedling Road, Washington,
England, NE38 9BZ
www.rocketmedical.com

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

This device is not
manufactured
with natural
rubber latex

Rx ONLY

Potential of DEHP Exposure: Medical procedures using PVC medical devices containing DEHP have the potential to lead to DEHP exposure due to DEHP leaching from the device. The extent of exposure largely depends upon the medical treatments administered, the duration of the treatment, and in the case of plastic blood bags, by the length of storage and the storage temperature. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilize mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure. (Full statement at www.rocketmedical.com).

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State. Please contact pncf@rocketmedical.com

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Copyright© 2013-2022