# Rocket FBS Lactate Kits

## INSTRUCTIONS FOR USE

**Indications:** These instructions apply to Rocket Fetal Blood Sampling Kits for Lactate measurement, product code R57028-00-SW. Procedure kit for single patient use, to obtain up to 5 fetal blood lactate samples in cases of suspected fetal distress during labour. If additional sampling is necessary or if labour is protracted, a further pack may be required. This device should only be used by or under the supervision of trained personnel and in conjunction with clinical practice guidelines.

**This kit is indicated for use with handheld lactate analysers that are compatible with the collection of blood in capillary tubes containing Na/Li unbalanced heparin in high concentration (240USP/ml).**

**Contraindications:** This procedure is contraindicated with placenta previa, when identification of the presenting part is uncertain, in the presence of genital infections (e.g. herpes, Group B streptococcus, gonorrhoea, HIV/AIDS and Hepatitis B) or where the mother is a confirmed carrier of haemophilia and the fetus is either affected or of unknown status, or acute fetal compromise with associated fetal bradycardia of >3 mins.

**WARNING: Do not use this kit without prior verification that it is suitable for use with the analyser.**

**CAUTION: Avoid fetal face, fontanels and genitalia.**


**CAUTION: DO NOT USE ETHYL CHLORIDE SPRAY – this may adversely affect the internal surface of the plastic amnioscope.**

**CAUTION: Do NOT use the sample if the blood clots**

**CAUTION: Do not shake the filled collection tube or use wires and magnets to mix heparin. Use of these techniques may cause air entrainment into the sample. Samples containing air bubbles may be rejected by the analyser.**

### Preparation:

- Using aseptic technique, remove the kit from outer bag and unfold outer drape to form a sterile field.
- Follow local hospital procedure to prepare the perineum and vagina using antibacterial solution.
- Select a capillary tube from the kit and load into the holder. Locate the end of the tube into the holder at ❶ and then gently advance into the rear tube grip at ❷. Ensure the tube is secure before sampling.
- The sampling wand comes pre-fitted with an FBS blade. Lightly holding the shaft and rear finger grip (A), push the grip forward to expose the blade (B). **Ensure the blade is locked forward before use.**
- Locate the Amnillum over the slot in the endoscope and push gently to engage. The light will activate when the light source is correctly fitted. Remove by pushing off with the thumb.
- Insert endoscope complete with obturator into the vagina against the fetal presentation. Remove obturator. Ensure that the cervix is not trapped. Position the bevelled end by rotating endoscope to obtain the best possible seal against the fetal presentation.
- Clean the puncture area with the x-ray detectable gauze swabs and forceps provided.
- With a clean swab stick, apply a *thin* layer of petroleum jelly (supplied) to the puncture area to promote blood aggregation.
- IMPORTANT:** Remove excess petroleum jelly with a new dry swab stick before skin puncture.
- Perform skin puncture through the endoscope using the FBS blade holder or Sampling Wand. Either ❶ make a vertical stab, rotate the blade 90° and make a second stab to form a cross or ❷ make a vertical stab, then angle the blade upwards and extend using the front edge of the blade to produce a single incision 4-5mm long. Always wait for a good sized (3-4mm) bead to form on the surface before collection. Never slash at the scalp and avoid multiple small stabs. Rocket® FBS blades have a maximum depth of penetration of 2.2mm.
- Collect the blood sample using the capillary tube. Typically, handheld lactate analysers require up to 10µl of blood, which can be obtained by filling the distal 1-2cm of the sampling tube only.
- Attach a Rocket Capillary Tube Dispenser to the proximal end of the tube.
- Following the analyser manufacturer's instructions, insert a test strip into the analyser and squeeze the bulb to dispense a suitable bead of blood onto the test strip. Complete the analysis.

### Post sampling:

- Haemostasis is accomplished by applying pressure to fetal puncture site with a swab for one full contraction.
- If test results and local clinical guidance so indicate, repeat procedure.

**For Single Use Only.** Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross contamination and infection with transmissible diseases.

### Disposal:

This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

Rx ONLY

STERILE EO

ROCKET MEDICAL PLC Sedling Road, Washington, England, NE38 9BZ [www.rocketmedical.com](http://www.rocketmedical.com)

EC REP

Rocket Medical GmbH  
Am Rosengarten 48,  
15566 Schöneiche.  
Germany

**DO NOT RESTERILISE**

Unless opened or damaged, contents of package are sterile.

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Copyright©2013-2020
All Rights Reserved (GB)