
## Instructions for Use


**Scope:** These instructions cover all Rocket Pro-Lux Sigmoidoscopes and derivatives.

**Indications:** This device should only be used by or under the supervision of trained personnel and in conjunction with current local clinical practice guidelines.


**Disposal:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment. Whilst not required by the manufacturers, where national, local or hospital trust policies specifically require the removal of battery cells prior to disposal, the instructions below should be followed.

1. Remove from packaging
2. Check that the light illuminates before patient use. Press the button on the Pro-Lux handle to operate the lightsource
3. Attach bellows tube onto connector as shown.
4. Lubricate well as required
5. Remove obturator following insertion

6. Press the button on the Pro-Lux handle to operate lightsource
7. Ensure lens cover is tightly closed
8. Following procedure, remove grey cap to expose battery pack ready for removal prior to disposal.

**HAZARD: Battery packs must NOT be allowed to come into contact with the batteries of other devices as they may short-out, leak fluid, overheat or explode.**


9. To ensure safety, place each battery pack into a separate cell of the disposal tray provided with each carton of Pro-Lux Proctoscopes
10. Dispose of the filled tray in accordance with local hospital policy.

This device is not manufactured with natural rubber latex

ROCKET MEDICAL PLC Sedling Road, Washington, England, NE38 9BZ [www.rocketmedical.com](http://www.rocketmedical.com)

Rocket Medical GmbH  
Am Rosengarten 48,  
15566 Schöneiche,  
Germany

### NON STERILE

Manufactured in the UK by Rocket Medical plc. 2-4 Sedling Road, WASHINGTON, Tyne & Wear, NE38 9BZ.

#### Potential of DEHP Exposure

Medical procedures using PVC medical devices containing DEHP have the potential to lead to DEHP exposure due to DEHP leaching from the device. The extent of exposure largely depends upon the medical treatments administered, the duration of the treatment, and in the case of plastic blood bags, by the length of storage and the storage temperature. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure. Full statement available [www.rocketmedical.com](http://www.rocketmedical.com)

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


ROCKET MEDICAL PLC Sedling Road, Washington, England, NE38 9BZ [www.rocketmedical.com](http://www.rocketmedical.com)

**Rocket Medical GmbH**  
Am Rosengarten 48,  
15566 Schöneiche.  
Tyskland

### USTERIL

Fremstillet i UK af Rocket Medical Plc. 2-4 Sedling Road, WASHINGTON, Tyne & Wear, NE38 9BZ.

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