INSTRUCTIONS FOR USE

ZDOCK049 Rev. 19

	Page
UK	2-3
DA	4-5
DE	6-7
ES	8-9
FI	10-11
FR	12-13
HU	14-15
IT	16-17
NL	18-19
PT	20-21
RO	22-23
AR	24-25
RU	26-27
SV	28-29

HSG and HyCoSy Catheter INSTRUCTIONS FOR USE

Catheter Product Codes: R57422, R57422-00-ST, R57424, R57424-00-ST

Procedure Pack Product Codes: R57422-00-PK, R57422-ST-PK, R57424-00-PK, R57424-ST-PK

Intended Use: For the instillation of aqueous x-ray and ultrasound contrast media during HSG and HyCoSy procedures.

Indications: To aid in the diagnosis of fallopian tube blockage, uterine abnormality and other pathological conditions associated with primary and secondary infertility. This device is only for use by or under the supervision of trained personnel and in accordance with local and national guidelines.

Contraindications: Not for use in patients with a vaginal infection, sexually transmitted disease, active vaginal bleeding or suspected or confirmed pregnancy.

WARNING: Maximum balloon volume = 2.0 ml.

Do not over inflate the balloon as this may cause severe discomfort, or balloon rupture and loss of the contrast seal.

Only suitable for use with aqueous contrast media. Do not use with oil-based contrast as this will cause the balloon to rupture.

CAUTION: NOT for use with pressure injectors.

CAUTION: The balloon should be protected from light until use.

Preparation:

1. Use aseptic technique throughout. Follow local hospital policy regarding PPE, which should include provision for the safe handling of contrast media. This diagnostic procedure can normally be carried out without medication, although light IV sedation and/or analgesia may be necessary in very anxious patients.
2. Open the package and remove the device and syringe. If not pre-attached, connect the syringe to the shorter side arm of the catheter (balloon port). Inject 2.0ml of air into the balloon and lightly palpate with the fingers to confirm balloon patency. Withdraw **all** residual air from the balloon.
3. Place the patient into the lithotomy position with the knees raised, expose the cervix with a Cusco or Sims type speculum. Assess the orientation of the cervix and establish the position of the external cervical os. Swab away any excess mucus with sterile normal saline or an aqueous antiseptic solution.
4. Pre-form the catheter with the wire stylet in place (applies only to R57422-00-ST, R57424-00-ST, R57422-ST-PK and R57424-ST-PK versions with integral forming wire).
5. Pre-fill the catheter with contrast media to avoid the risk of erroneous diagnosis caused by air bubbles.
6. To aid insertion, a small quantity of water-soluble lubricant may be applied to the catheter tip and balloon.
7. Place the distal tip of the catheter into the external cervical os. Using a very gentle rotating motion, advance the catheter into the uterus, avoiding the uterine wall. The patient may report slight discomfort when the catheter passes through the internal os and into the uterine cavity. The outer sliding sheath may be used to provide extra rigidity to the catheter to aid initial insertion.
8. Using the supplied syringe, inflate the balloon with up to **2.0ml** of air, sufficient to create a seal (Figure. 1). Gently withdraw the catheter until it is located against the internal os. The pressure of the balloon against in the internal os creates a seal and prevents leakage of the contrast media.
9. The catheter may also be inflated in the endocervical canal if insertion beyond the internal os is deemed undesirable.
10. Attach a 10 or 20ml Luer Lock syringe to the rear port, ensuring the tap is in-line with the catheter. Slowly instil 10-20ml of contrast media, sufficient to adequately demonstrate the uterine cavity and fallopian tubes. Further volumes may be added by turning the tap through 90° to the 'Off' position and replenishing the syringe.


Figure 1

To remove the device: Attach the supplied syringe to the balloon port, open the tap and withdraw the plunger to deflate the balloon. Gently withdraw the catheter from the uterus.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.


This device is not manufactured with natural rubber latex

Rx ONLY

STERILE

ROCKET MEDICAL PLC
Sedling Road, Washington,
England, NE38 9BZ
www.rocketmedical.com

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

DO NOT RESTERILISE

For Single Use Only. Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases. **Unless opened or damaged, contents of package are sterile.**

Distributed in the USA by **ROCKET MEDICAL**

50 Corporate Park Drive #890, PEMBROKE, MA 02359. USA Tel: +1 781 749 6223

For HSG Procedure Packs Only (Product Codes R57422-ST-PK, R57422-00-PK, R57424-ST-PK and R57424-00-PK:

Contents:

- Catheter (Rocket Medical)
- Filling Kwill (Rocket Medical)
- Uterine Sound (Rocket Medical)
- Luer Lock Extension Tube (Rocket Medical)
- Two-piece Adhesive Aperture Drape with Fenestration (Rocket Medical)
- Abdominal Drape (Rocket Medical)
- Gauzes x 5 (Bunzl Healthcare)
- Dilator (Marina Medical)
- Optilube Lubricant (Optimum Medical)
- Speculum (Pelican Feminine Healthcare)
- 20mL Syringe (Becton Dickinson)
- 2.5mL Syringe (Terumo)
- Spongeholding Forceps (Pelican Feminine Healthcare)
- 4 x Swabsticks (Applimed SA)
- Blue Applicator (Promepla SA)
- Small Pot for Skin Prep (not to be used for contrast media or medicinal agent) (Cups Direct)

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