## Ultra-Fine™ Diathermy Electrodes

Product Codes: R57800, R57801, R57802, R57803, R57804, R57807, R57851, R57852 and R57860

### INSTRUCTIONS FOR USE

#### Scope:

These instructions cover all Ultra-Fine™ LLETZ Loops, Ultra-Fine™ DELTA™ LLETZ Loops, Ultra-Fine™ Coagulation Balls and Ultra-Fine™ NETZ Needles. (Figure 7).

#### Indications:

For large loop excision or needle excision of the transformation zone in the diagnosis and treatment of Cervical Intraepithelial Neoplasia (CIN i,ii,iii) or cervical dysplasia. These devices should only be used by or under the supervision of personnel with recognised training in the treatment and management of CIN such as that approved by the BSCCP and the Royal College of Obstetricians & Gynaecologists.

#### Contraindications:

These devices are contraindicated for use in patients with suspected or confirmed pregnancy, invasive cancer that is visible on examination, acute or active inflammation or infection of the cervix, endometrium, fallopian tube, ovary or peritoneum (cervicitis, endometritis, tubo-ovarian inflammatory disease or pelvic inflammatory disease), known or suspected cervical changes secondary to diethylstilboestrol (des) intrauterine exposure, patients with bleeding disorders or patients less than 3 months post-partum.

#### WARNINGS:

Do not allow the wire loop to come into contact with metal speculae or other metal instruments. Contact while activated will cause the wire loop to burn out. Use of higher than recommended power settings will cause the wire loop to burn out.  
 Do not use in the presence of flammable gases or liquids, or objects, such as cotton or gauze in oxygen-enriched environments. Do not use in the presence of flammable gases that have been released from the bowel. These may be ignited by sparks produced during the normal use of HF surgical equipment.  
 There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Avoid pooling of flammable liquids  
 Do not use in the presence of nitrous oxide or other oxidizing agents.  
 The electrode will become hot during use – keep away from flammable materials.  
 For patients with electrically conductive implants, a possible hazard exists due to concentration or redirection of high frequency currents, and when in doubt, qualified advice should be obtained.  
 The electrode must be fully engaged with the hand switch, ensuring that no bare metal remains exposed. If not fully engaged, and the user makes contact with exposed bare metal, there is a risk of electric shock when the electrode is activated.  
 If a poor electrical connection is made, discard the electrode. Select a replacement electrode to complete the procedure.  
 Do not use if the pack is opened or damaged – discard and complete the procedure using another device.  
 This IFU should be read in conjunction with the user manual for the electrosurgical generator in use. It is not intended as a replacement for comprehensive professional training.

#### Recommended Operative Technique:

(1) Following colposcopic examination and establishing the presence of CIN or associated cervical conditions, infiltrate the cervix with a local anaesthetic (such as 2% lignocaine solution BP) using a dental syringe. (Figures 1 and 2).

(2) Insert the selected diathermy loop or needle into an approved finger switch/probe and connect to a surgical diathermy unit following the manufacturer's recommendations. The selection of the electrode size and shape is dependent upon:

- (a) Lesion size and type
- (b) Depth of excision
- (c) Personal choice and surgical technique.

(3) Using Table 1 as a guide, set the power output following the manufacturers recommended operating instructions.

(4) LLETZ: Pass the loop horizontally into the cervical stroma to the desired depth and then in a downward sweep complete the excision. If the lesion is very large then two or more sweeps may be made on either side of the os until the process has removed the entire diseased area (Figures 3 and 4).

(5) NETZ: Insert the needle to the desired depth. Pass the needle through the tissue to cut in a manner similar to that of a scalpel.

**IMPORTANT:** The wire of the loop or needle does not play any part in the 'cutting' of the tissue. The incision is created *solely* by the passage of the high frequency current around the wire. Deformation of the wire indicates forced passage through tissue.

(6) The excised tissue or 'biopsy' should be sent for histological examination to ensure that the whole abnormality has been eradicated (Figure 5).

(7) Swab the resulting crater of any free blood and apply a ball electrode to coagulate any bleeding points. Ensure that the crater is left dry (Figure 6).

(8) Following hospital or clinical preferences, the patient may be advised to refrain from intercourse for three weeks and to expect a light bloody discharge for 4 - 5 days. An antiseptic vaginal cream may be provided to help prevent infection and to soothe any soreness.

(9) The patient should be returned for follow up colposcopy at 6 or 12-months' post treatment for 'test of cure' repeat cytology, or as per local guidelines.

**Disposal:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, 'those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

Rx ONLY

STERILE EO

**ROCKET MEDICAL PLC**  
 Sedling Road, Washington,  
 England, NE38 9BZ  
[www.rocketmedical.com](http://www.rocketmedical.com)

EC REP

**Rocket Medical GmbH**  
 Am Rosengarten 48,  
 15566 Schöneiche.  
 Germany

#### DO NOT RESTERILISE

Unless opened or damaged, contents of package are sterile. Sterile unless package has been opened or damaged.  
 Classification: BF.

Ultra-Fine™, DELTA™ and NETZ™ are trademarks of **ROCKET MEDICAL PLC**.

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#### Associated Devices

The diathermy electrodes are intended to be used with a diathermy machine capable of the following outputs:

Modes: cut, coagulation, fulguration.  
Output voltage: 3200V max.  
Output Power: 0-60W  
Output Frequency: 500KHz (monopolar mode)

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(AR)

## Recommended Diathermy Output Settings (Table 1)

(UK)

Loop Size	Blend	Watts	Example
Small/Square	BLEND 1	25-30	Surface work on cervix
Small/Square	BLEND 1	25-30	TZ excision
Medium including DELTA™	BLEND 1	30-40	TZ excision
Large including DELTA™	BLEND 2 or 3	30-40	TZ excision
3mm Ball	COAG	35-40	Crater coagulation
5mm Ball	COAG	35-55	Crater coagulation
NETZ needles	BLEND 1	25-30	Needle excision

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Rev.16 Copyright© 1999-2020
All Rights Reserved (GB)
Copyright© 1999-2020