
# Embryon<sup>®</sup> BulbTip<sup>™</sup> Embryo Transfer Catheters

## INSTRUCTIONS for USE

### 1. SCOPE:

These instructions apply to Rocket Medical BulbTip<sup>™</sup> Embryo Transfer Catheters and Trial Catheters, Product Codes R57635-00-18, R57635-00-23, R57635-EC-18, R57635-EC-23, R57636-00-18 and R57636-00-23.

### 2. INDICATIONS:

For the placement of human embryos into the uterine cavity. This device may only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines.

### 3. CONTRAINDICATIONS:

Not for use in patients with a vaginal infection, a sexually transmitted disease or suspected pregnancy.

#### **WARNING: CATHETERS WITH ECHOGENIC BANDS:**

**Do NOT withdraw the inner catheter through the outer sheath while the tip of the catheter is in the uterus or endocervical canal. Detachment of the echogenic band may result.**

**Remove both inner and outer catheter together and recommence the insertion procedure**

### 4. DESCRIPTION:

The device is constructed of two primary components

1. Soft inner catheter
2. Outer sheath with 5 x 1cm markings from the tip to indicate depth of penetration into the endocervix

The transparent finger grip has 5 ribs moulded at 1cm intervals. When the inner catheter is inserted into the outer sheath, the coloured marker bead will indicate the position of the catheter tip in relation to the end of the outer sheath. Position of the coloured bead corresponds with:

- a) Proximal edge - tip protrusion = 0cm (flush with the end of sheath)
- b) 1st proximal rib - tip protrusion = 1cm
- c) 2nd proximal rib - tip protrusion = 2cm.
- d) Catheter fully inserted - tip protrusion = 4cm

BulbTip<sup>™</sup> Echogenic Transfer catheters have an echogenic band applied 2mm from the distal tip. This echogenic marker may aid placement with ultrasound. Accidental detachment of the echogenic band from the inner has been reported in <0.016% of cases.

### 5. INSTRUCTIONS FOR USE:

**Cervical Assessment: Rocket Medical plc.** strongly recommends that the position and patency of the cervix and canal is assessed 2-3 weeks before the date of initial transfer. If more than 3 months has elapsed since the previous transfer or assessment, then a further examination should be considered. The total uterine length, including endocervix, should be established by uterine sounding or ultrasonography at the initial patient consultation.

Place the patient in the lithotomy position and expose the cervix using a conventional Cusco or Sims speculum. The procedure can normally be carried out without medication, although light sedation may be necessary in very anxious patients.

Assess the position of the cervix and establish the position of the external cervical os. Swab away any excess mucus with normal saline solution. Using aseptic technique, prepare a matching **Trial Transfer Catheter** and follow one of the two techniques below:

- a) **No significant anteversion or retroversion, cervical os fully visible:** Assemble the inner trial catheter into the outer sheath and secure the luer fitting into the clear outer handle. Locate the distal tip of the trial catheter into the external cervical os and with a very gentle rotating motion, advance the catheter into the uterus. The patient may report slight discomfort when the catheter passes through the internal os and into the endometrial cavity. If the trial transfer catheter will not pass freely into and through the endocervical canal, follow the technique below.
- b) **Anteversion, retroversion with cervical os not fully visible or previous catheter insertion difficult or obstructed:** Remove the inner trial catheter from the outer sheath. Gently grip the outer sheath between gloved forefinger and thumb and *gently* increase the curve to match the amount of anteversion or retroversion. Do not bend sharply as this will kink the outer sheath.

Identify the position of cervical os and advance the outer sheath 2 x 1cm markings. Insert the inner catheter into the clear handle and advance until the inner catheter passes up to the internal os.

Observing the relative position of the marker bead and the moulded finger grooves will indicate the catheter tip position, from the distal end of the outer sheath.

Patients who have a fully obstructed cervix or impassable endocervical canal should be considered for cervical surgery or alternative fertility treatment before further attempts at endocervical transfer are made.

**Drawing up embryos:** Select a transfer catheter to suit the needs of the patient and surgical preference.

- a) Using aseptic technique, remove the device from the packaging and attach a sterile microsyringe. Do not remove the inner catheter from the outer sheath as this provides additional support and control.
- b) Under microscopic control, draw 100µl of prepared medium to fill the catheter. Draw up a further 30µl containing the embryos for transfer. The embryos should be retained in the distal third of the catheter until the moment of injection into the uterus. Remove the microsyringe, expel any residual air and fill with medium ensuring that no air remains in the syringe.

OR:

- c) Under microscopic control, draw 50µl of prepared medium into the distal third of the catheter. Draw up 10-20µl of air. Draw up a further 30µl containing the embryos for transfer. The embryos should be retained in the distal third of the catheter until the moment of injection into the uterus. Remove the microsyringe and fill with medium, ensuring that no air remains in the syringe.
- d) Maintaining aseptic technique, transfer the complete assembly to the operating room immediately prior to transfer.

**Transferring the embryos:**

- a) Ensure that trial transfer and cervical assessment have been completed before attempting a 'live' transfer. Place the patient in the lithotomy position and expose the cervix using a conventional Cusco or Sims speculum. The procedure can normally be carried out without medication, although light sedation may be necessary in very anxious patients.
- b) Transvaginal or transabdominal ultrasound may be useful in establishing correct placement of the catheter tip in the uterus.
- c) Assess the position of the cervix and establish the position of the external cervical os. Swab away any excess mucus with normal saline solution or prepared medium.
- d) Follow the techniques described for the placement of the trial Transfer Catheter to place the loaded embryo transfer catheter.
- e) If a trial transfer is completed immediately before live transfer, the outer sheath of the trial catheter should be left in position and the trial inner removed. The inner loaded embryos should be substituted and inserted to the measured depth indicated by the marker beads.
- f) Once in position, inject approximately 20µl of medium to expel the embryos.
- g) Whilst maintaining pressure on the syringe plunger, gently withdraw the catheter from the uterus and return it to the laboratory to check microscopically that all the embryos have been expelled.
- h) The patient may then be returned to the ward for subsequent discharge.

**DISPOSAL:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

**STERILE**

**Rx ONLY**

**MD**

**UK  
CA**

**ROCKET MEDICAL PLC**  
Sedling Road, Washington,  
England, NE38 9BZ  
[www.rocketmedical.com](http://www.rocketmedical.com)

**EC REP**

**Rocket Medical GmbH**  
Am Rosengarten 48,  
15566 Schöneiche.  
Germany

**Unless opened or damaged, contents of package are sterile**


Distributed in the USA by Rocket Medical. 50 Corporate Park Drive #890, PEMBROKE. MA 02359. USA Tel: +1 781 749 6223

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ZDOCK67 2022-02-18 Rev.24 Copyright© 2009-2022. **ROCKET MEDICAL PLC** All Rights Reserved. (GB)
ZDOCK67 2022-02-18 Rev.24 Copyright© 2009-2022. **ROCKET MEDICAL PLC** Alle Rechte vorbehalten. (DE)
ZDOCK67 2022-02-18 Rev.24 Copyright© 2009-2022. **ROCKET MEDICAL PLC**
ZDOCK67 2022-02-18 Rev.24 Copyright© 2009-2022. **ROCKET MEDICAL PLC** Tous droits réservés. (FR)
ZDOCK67 2022-02-18 Rev.24 Copyright© 2009-2022. **ROCKET MEDICAL PLC** Tutti i diritti riservati. (IT)
ZDOCK67 2021-02-18 Rev.24 Copyright© 2009-2022. ROCKET MEDICAL PLC