Rocket® Single and Double Lumen Oocyte Aspiration Needles INSTRUCTIONS FOR USE

Scope: These instructions cover all Rocket® Single and Double Lumen Oocyte Aspiration Needles and derivatives

Indications: For the aspiration and flushing of ovarian follicles under real time ultrasound visualisation. This device should only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines such as those published by the National Institute of Clinical Excellence and the Human Fertilisation & Embryology Authority.

Contraindications: These devices are contraindicated for use in patients with vaginal infections, sexually transmitted diseases, suspected pregnancy and where the pelvic anatomy cannot be accurately assessed by ultrasonography.

Preparation: **ROCKET MEDICAL PLC** strongly recommends confirmation of the position of the uterus and pelvic organs by ultrasonography prior to the procedure. Similarly, the presence of developed follicles should be confirmed by serial ultrasound during the ten days prior to the procedure. The procedure is normally carried out under local analgesia and/or supplemented light sedation. General anaesthesia may be necessary in very anxious patients.

Follicular Aspiration: It is essential for a successful outcome that oocytes are harvested with the least possible trauma to both patient and ova.

Using aseptic technique:

- Following the manufacturers guidelines, cover the transvaginal ultrasound probe with a sterile sheath. A small amount of ultrasound gel may be added to the sheath to improve picture definition
- Select a sterile needle guide approved for use with the vaginal ultrasound probe and assemble following the manufacturers instructions.
- Open the needle pack carefully, taking care not to touch or damage the needle bevel. The use of damaged needles will cause increased discomfort to the patient and may result in loss of the oocyte.
- For double lumen needles/flushing technique: prior to inserting the needle into the needle guide, attach a syringe of flushing media to the flushing port, open the tap and flush the channel with 2ml of media. Close the 2-way tap. For single lumen needles: prime (if required) using the connection port on the underside of the bung/stopper.
- Attach a sterile sample (Falcon) tube to the bung/stopper.
- Connect the tubing connector to a low power vacuum pump such as the R29.655 Rocket® Craft™ Pump.
- For double lumen needles/flushing technique: set the pump to deliver a vacuum of 200mmHg. Aspirate 2-5ml of flushing media into the Falcon tube and discard.
- Place the needle into the needle guide following the manufacturer's instructions. Figure 1 shows a common probe and needle assembly.
- Set the aspiration vacuum level to your personal preference using the table as a guide.
- Ensure all connections are air tight and the tube set is free from any constrictions as these cause turbulences which greatly increases the risk of oocyte damage.
- With the patient in lithotomy, introduce the needle/probe assembly into the vagina, advancing into the anterior fornix to visualise each ovary. Once the ovaries have been identified, introduce the needle and advance through the vaginal wall and ovarian stroma into the target follicle.
- Activate the vacuum pump and aspirate the follicular fluid into the sampling tube. Check microscopically for the presence of an oocyte. For double-lumen needles, if desired, 2-5ml of media can be injected via the flushing limb (after opening the 2-way tap) to distend the follicle prior to subsequent aspiration. When using single lumen needles, follicular flushing can be performed via the flushing port on the underside of the bung/stopper. Complete the harvesting from one ovary before commencing on the other.
- If the needle becomes blocked during the procedure, rotate the needle within the follicle to ensure that the needle is not blocked by follicular wall tissue

Single lumen blockage: if the needle remains blocked, remove the needle from the patient and flush through with media using the flushing port on the underside of the bung/stopper.

Double lumen blockage: the needle can be flushed without removal from the ovary.

Check that the needle is clear by aspirating media through the needle set. If the blockage is still present, high vacuum aspiration (400-500mmHg) may be performed with the needle outside the patient. **WARNING:** the use of high suction when the needle is inside the ovary should not be performed. This may lead to turbulent flow, which may damage the oocyte.
- Flushing may be performed following harvesting to clear the needle prior to oocyte identification.


Figure. 1

Tube Set Length	Single Lumen		Double Lumen	
	16G	17G	16G	17G
55cm	80	110	130	150
70cm	90	130	150	170
90cm	100	150	170	190

This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, 'those pertaining to human health & safety and care of the environment

This device is not manufactured with natural rubber latex

Rx ONLY

STERILE EO

ROCKET MEDICAL PLC
Sedling Road, Washington,
England, NE38 9BZ
www.rocketmedical.com

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

For Single Use Only. Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases.

Unless opened or damaged, contents of package are sterile.

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Copyright© 2011-21.
All rights reserved.
ZDOCK122 Rev.17 2021-01-18 Copyright© 2011-21. **ROCKET MEDICAL PLC** Todos los derechos reservados. (ES)
ZDOCK122 Rev.17 2021-01-18 Copyright© 2011-21. **ROCKET MEDICAL PLC** Tutti i diritti riservati. (IT)
ZDOCK122 Rev.17 2021-01-18 Copyright© 2011-18. **ROCKET MEDICAL PLC** Todos os direitos reservados. (PT)