Seldinger Chest Drain Kit

INSTRUCTIONS FOR USE

Scope: These instructions cover all R54544-12-00 and R54544-18-00 Seldinger Chest Drain Kits.

The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques. It is recommended that British Thoracic Society Guidelines for the insertion of chest drainage tubes are followed.

Upon removal from package, inspect the product to ensure no damage has occurred.


Product Description: The Rocket® Seldinger Chest Drain Kit contains a 10ml syringe, introducer needle, wire guide, 10cm dilator(s), fir tree and wide bore tubing set connectors, 3 way tap, coaxial chest tube inserter and clear drainage catheter with radiopaque stripe.

Indications: The product has been designed for percutaneous introduction of a chest tube for drainage of air and/or pleural fluid.

Contraindications: For drainage of pleural space only
12fg – not for drainage of thick fluids.

Instructions For Use:

1. Following local hospital policy, prepare the chest tube insertion site with an approved solution such as povidone-iodine and drape to maintain and aseptic technique.
2. Identify and mark the insertion point following BTS Guidelines, with specific attention paid to the 'safe-triangle'.
3. After injection of local anaesthetic, make a small skin incision 4-5mm.


CAUTION: Insertion should be just ABOVE the rib to avoid damage to the intercostals neurovascular bundle.

The introducer needle has changed from a Tuohy needle and is now a straight Seldinger needle as below. When removing the protective sheath, hold the silicone ring to prevent it from being removed from the needle.


The measurement provided by using the silicone ring can be used to set the position of the sliding stop on each dilator to suit the size of the patient and help prevent over insertion of the dilators.

4. Attach the introducer needle to a syringe filled with a local anaesthetic e.g. lignocaine, and administer additional local anaesthetic to the skin and underlying tissues.
5. Advance the introducer needle over the *superior* border of the rib and into the pleural space. Fluid or air should be aspirated to verify correct position. The needle should be introduced and directed with appropriate orientation inferiorly or superiorly.
6. Remove the syringe and advance the soft "J" end of the wire guide through the needle and into the pleural space. Before removing the needle, move the silicone measurement ring down to the patient's skin.

CAUTION: The wire guide should pass through the needle and advance into the pleural space without resistance.


7. Remove the introducer needle.
8. Using the measurement obtained using the depth marker on the introducer needle, set the position of the sliding stop on each dilator to suit the size of the patient.


WARNING: Over insertion of the dilator into the chest cavity is associated with serious injury and **MUST** be avoided.
Only insert the dilator sufficient to dilate the superficial track.

CAUTION: Ensure sufficient length of wire guide remains external to the chest to facilitate controlled introduction of the dilators and drain.

9. While maintaining the wire guide position, dilate the track and opening into the pleural space by advancing and withdrawing, one by one, the supplied dilator(s) over the wire guide. Introduction into the pleural space is facilitated by rotating and advancing the dilators in the same plane as the wire guide to prevent kinking.
10. Remove the final dilator from the wire guide


11. Locate the drain set. Fully insert the stiffening tube, tighten the valve if required to ensure the stiffening tube remains fully inserted to the distal plug. Insert the free end of the guidewire into the distal catheter plug ensuring it travels through the stiffening tube.

12. With the wire guide still positioned within the pleural space, advance the chest tube inserter/chest tube assembly over the wire guide and into the pleural space.

CAUTION: It is important to advance the chest tube assembly into the pleural space whilst maintaining its direction in the same plane as the wire guide. This will make introduction easier and avoid kinking of the wire guide.

CAUTION: Ensure all side holes of the chest tube are positioned within the pleural space.

13. Remove the wire guide and chest tube inserter leaving the chest tube in place.
14. The chest tube can now be sutured to the skin and is ready for connection to a suitable sealed drainage system such as the Rocket® Blue Bottle or Rocket® Ambulatory Bag
15. Confirm the correct positioning of the drain with imaging as soon as possible.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

Rx ONLY

STERILE EO

PHT
DEHP

ROCKET MEDICAL PLC Sedling Road, Washington,
England, NE38 9BZ www.rocketmedical.com

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

CONTINUOUS USE SHOULD NOT EXCEED 28 DAYS.

For Single Use Only. Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases. **Unless opened or damaged, contents of package are sterile**

Potential of DEHP Exposure

Medical procedures using PVC medical devices containing DEHP have the potential to lead to DEHP exposure due to DEHP leaching from the device. The extent of exposure largely depends upon the medical treatments administered, the duration of the treatment, and in the case of plastic blood bags, by the length of storage and the storage temperature. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure.

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