INSTRUCTIONS FOR USE

Water Trap Set for CRAFT™ Suction Pumps

Scope: These instructions cover: R57685 Water Trap Set for CRAFT™ Suction Pumps.

Indications: For the connection of oocyte aspiration needles to the CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps permitting vacuums up to 400mmHg⁻¹ to be applied for the recovery of oocytes as part of the treatment of infertility. This device should only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines.

Contraindications: Do not allow set to become contaminated with fluid. If set becomes contaminated replace immediately.

Set Up Procedure:

1. Remove and discard reusable bottle set if present, including bottle, grey stopper, metal tubes, silicone tubing and adapter (applies to older models only).
2. Observing aseptic technique, open the Water Trap Set for CRAFT™ Suction Pumps onto the sterile field and pass the bottle out to the circulating assistant. Retain the distal end of the long patient connection tube in the sterile field.
3. Carefully insert the bottle into the pump clip.
4. Locate the short length of tube and attach firmly to the vacuum outlet indicated on the diagram. Ensure there are no kinks or restrictions in the tube run.
5. Fully occlude the *patient* connection tube between forefinger and thumb downstream of the water trap bottle.
6. Ensure the pump is connected to the correct power source.
7. Maintain the occlusion and operate the Standard Vacuum footswitch.
8. Rotate the vacuum control to set the vacuum to between 100-200mmHg.
Note: Collection needles with tube sets >70cm in length typically require vacuums >150mmHg to achieve satisfactory flow rates.
9. Once set, operate the HIGH vacuum footswitch (or occlude the red button) to achieve maximum vacuum. The vacuum reading should exceed 400mmHg⁻¹.
10. Release the tube set, the system is now ready for use.
11. Do NOT allow the bottle to fill with fluid. Replace immediately if the bottle becomes contaminated.

CAUTION: Vacuum settings should always be made with the patient connection tube fully occluded.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

STERILE EO

PHT
DEHP

Rx
ONLY

ROCKET MEDICAL PLC
Sedling Road, Washington,
England, NE38 9BZ
www.rocketmedical.com

MD

UK
CA

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

For Single Use Only: Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases. **Sterile unless package has been opened or damaged. NOT FOR LONG TERM USAGE.**

Potential of DEHP Exposure: Medical procedures using PVC medical devices containing DEHP have the potential to lead to DEHP exposure due to DEHP leaching from the device. The extent of exposure largely depends upon the medical treatments administered, the duration of the treatment, and in the case of plastic blood bags, by the length of storage and the storage temperature. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure. (Full statement at www.rocketmedical.com).

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