## Rocket® Copeland® Fetal Scalp Electrode

**INDICATIONS:** A CTG electrode for attachment to the fetal presenting part to enable fetal heart rate monitoring. This device is only for use by or under the supervision of trained personnel and in accordance with clinical guidelines.

**CONTRAINDICATIONS:** Do not apply to the fetal face, fontanelles, suture lines or genitals. Not to be used in the presence of vaginal infection, HIV, infectious hepatitis, active herpes simplex virus. Not to be used in cases of known or suspected fetal coagulation disorder. Chronic fetal bleeding (e.g., due to partial placental abruption) leading to loss of fetal blood volume may result in a reduction in the margin of safety that the fetus can successfully respond to hypoxia.

### 1. PREPARATION

Conduct a vaginal examination to confirm that it is appropriate to use the fetal scalp electrode. Advance the Rocket® Copeland® between the examining fingers until the needle (A) is held pushed against the fetus by finger pressure on the serrated depression (B). Ideally apply to the scalp, avoiding the fontanelles and suture lines. The electrode may also be placed on the buttocks, legs or arms.

**DO NOT APPLY THE ELECTRODE TO THE FETAL FACE, THE FONTANELLES OR THE GENITALIA.**

### 2. APPLICATION AND CONNECTION

With the fingers of the other hand, rotate (C) approximately 180° ANTI-CLOCKWISE to open the needle and then approximately 180° CLOCKWISE to insert the needle into the fetal scalp. Do not rely on spring tension to insert. Firmly connect the fetal scalp electrode to the leg plate and the leg plate to the maternal attachment electrode. The leg plate cable will have a connector, which fits into the connecting parts of the fetal monitor. Follow the monitor manufacturer's instructions for methods of attachment.

### 3. REMOVAL

To remove, turn (C) ANTI-CLOCKWISE. It is not necessary to hold (A).

### 4. WARNINGS AND PRECAUTIONS

In the event of a poor trace or no trace at all, follow this sequence until a trace is established.

1. Check that you have followed the instructions and applied the fetal scalp electrode correctly.
2. Check that the fetal monitor and paper recorder is switched on.
3. Check that the maternal attachment electrode is correctly applied and is connected to the leg plate.
4. Check that the leg plate cable is correctly connected to the monitor.
5. Check that the fetal scalp electrode is firmly and securely inserted into the leg plate.
6. Re-position the electrode.
7. Replace the electrode.

**CAUTION Do not use if packaging is open or damaged**  
**CAUTION Do not over-rotate or apply excessive torque when applying or removing the device. This may cause damage to the spring/hook assembly, resulting in difficulty detaching the device.**

This device is not manufactured with natural rubber latex

**ROCKET MEDICAL PLC**  
 Sedling Road, Washington,  
 England, NE38 9BZ  
[www.rocketmedical.com](http://www.rocketmedical.com)


**Rocket Medical GmbH**  
 Am Rosengarten 48,  
 15566 Schöneiche.  
 Germany

### DO NOT RESTERILISE



For Single Use Only. Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases. Unless opened or damaged, contents of package are sterile

## Instructions for Use (GB) Brugsvejledning (DA) Bedienungsanleitung (DE) Istruzioni per l'uso (ES) Käyttöohjeet (FI) Mode d'emploi (FR) Instrucciones de uso (IT)

A) Gennemtrængende nåle  
 A) Eindringende Nadel  
 A) Aguja penetrante  
 A) Tunkeutuva neula  
 A) Ago di penetrazione  
 A) Aiguille de pénétration

B) Takkede fordybning  
 B) Genoppte Vertiefung  
 B) Depresión dentada  
 B) kohokuvioitua painaamaa  
 B) Avvallamento dentato  
 B) Dépression dentelée

ZDOCK167 2021-07-22 Rev.11 Copyright© 2021

**ROCKET MEDICAL PLC** All Rights Reserved **(GB)**

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ZDOCK167 2021-07-22 Rev.11 Copyright© 2021 (DA)
ZDOCK167 2021-07-22 Rév.II Copyright© 2021