Rocket Blunt Dissection Chest Drain Kit

INSTRUCTIONS FOR USE

Scope: These instructions cover all R54549-26-PK, R54549-28-PK and R54549-32-PK Rocket Blunt Dissection Chest Drain Kits and derivatives.

The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques. It is recommended that British Thoracic Society Guidelines for the insertion of chest drainage tubes are followed.

Upon removal from package, inspect the product to ensure no damage has occurred.


Product Description: The Rocket Blunt Dissection Chest Drain Pack contains a Disposable gown, fenestrated drape, skin preparation sponges and gallipot, forceps, No 11 Blade Scalpel, 10ml syringe, 20ml syringe, green needles, orange needle, tubing set connectors, clear drainage catheter with radiopaque stripe, suture and dressing pack.

Indications: The product has been designed for blunt dissection technique introduction of a chest tube for drainage of air and/or pleural fluid.


Contraindications: Non-pleural drainage.

Instructions For Use:

1. Following local hospital policy, and wearing the sterile gown included, prepare the chest tube insertion site with an approved solution such as povidone-iodine and drape with fenestrated drape included to maintain and aseptic technique.
2. Identify and mark the insertion point following BTS Guidelines, with specific attention paid to the 'safe-triangle'.
3. After injection of local anaesthetic e.g. lignocaine, make a small skin incision approx 10mm. and then continue administer additional local anaesthetic to the skin and underlying tissues
4. Allow the local anaesthetic to take effect. Use the scalpel to make an incision large enough that you will be able to place the catheter through it when you have created the path to the pleural space.

CAUTION: Insertion should be just ABOVE the rib to avoid damage to the intercostal neurovascular bundle.

5. By opening and closing the forceps make a path through the subcutaneous tissue and muscle into the pleural space.
6. Place the introducer in to the 4th fenestration on the catheter. Use the malleable introducer to support and position the drain as you insert the catheter into the chest. The catheter should be aimed apically for a pneumothorax, and basally for fluid
7. Insert the catheter gently through the chest wall, when the handle of the introducer nears the chest wall, remove it from the catheter so all the fenestrations will enter the chest.
8. Secure the drain by suturing the drain twice. The technique chosen should include adequate skin and subcutaneous tissue to ensure it is secure.

CAUTION: Ensure all side holes of the chest tube are positioned within the pleural space.

9. The chest tube can now be sutured to the skin and is ready for connection to a suitable sealed drainage system such as the Rocket Blue Bottle or Rocket Ambulatory Bag
10. Confirm the correct positioning of the drain with imaging as soon as possible.

How to use the Swann Morton Safety Scalpel

1. Grasp the scalpel and carefully extend the blade by moving the slide toward the tip of the scalpel, using the thumb of the hand holding the scalpel.
2. Extend the slider until you reach the positive stop, the slider will fit in the notch when it is completely extended.
3. To retract the blade, grasp the scalpel carefully and move the slide toward the back of the scalpel, using the hand holding the scalpel.
4. You should feel clicks as the blade is retracted and a positive stop once the blade is completely retracted.
5. To retract the blade permanently, move the slider past the notch at the back of the scalpel.
6. Dispose of the Swann Morton Retractable Disposable Scalpel in a puncture resistant container approved for sharps disposal, or in a manner consistent with your hospital procedure.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

Rx ONLY

STERILE EO

ROCKET MEDICAL PLC Sedling Road, Washington,
England, NE38 9BZ www.rocketmedical.com

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

CONTINUOUS USE SHOULD NOT EXCEED 28 DAYS.

Sterile unless package has been opened or damaged.

Store at room temperature. Avoid prolonged exposure to elevated temperature.

Potential of DEHP Exposure

Medical procedures using PVC medical devices containing DEHP have the potential to lead to DEHP exposure due to DEHP leaching from the device. The extent of exposure largely depends upon the medical treatments administered, the duration of the treatment, and in the case of plastic blood bags, by the length of storage and the storage temperature. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure.

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