# FBS Kits with Sampling Wand

## INSTRUCTIONS FOR USE

ZDOCKI95 Rev. 19

	Page
UK .....	2
DA .....	4
DE .....	6
ES .....	8
FI .....	10
FR .....	12
NL .....	14
NO .....	16
AR .....	18

# FBS Kits with Sampling Wand

## INSTRUCTIONS FOR USE

**Indications: Fetal Blood Sampling (FBS) Kit.** Procedure kit for single use, to obtain up to five fetal blood pH/gas samples in cases of suspected fetal distress during labour. If additional sampling is necessary or if labour is protracted, a further pack may be required. This device should only be used by or under the supervision of trained personnel and in conjunction with current clinical practice guidelines.

**Contraindications:** This procedure is contraindicated with placenta previa, when identification of the presenting part is uncertain, in the presence of genital infections (e.g. herpes, Group B streptococcus, gonorrhoea, HIV/AIDS and Hepatitis B), where the mother is a confirmed carrier of haemophilia and the fetus is either affected or of unknown status, prematurity (less than 34+0 weeks) or acute fetal compromise (for example, fetal bradycardia of >3 minutes).

**CAUTION: DO NOT USE ETHYL CHLORIDE SPRAY – this may adversely affect the internal surface of the plastic amnioscope.**

**CAUTION: Avoid fetal face, fontanelles or genitalia. Rocket Medical FBS blades have a maximum depth of penetration of 2.2mm**

**CAUTION: Failure to use correct blood volume, select correct analysis mode or the use of clotted samples or those contaminated with air bubbles may lead to errors in results or rejection of the whole sample. Carefully follow the analyser manufacturer's instructions.**



These instructions apply to all Rocket Medical Fetal Blood Sampling Kits with Sampling Wands. Contents of packs will vary. Rocket Medical supplies a range of FBS Kits suitable for a wide range of blood gas analysers. Please refer to the table below for the specification of the capillary tubes supplied within each kit and check your analyser manual and sampling volume to select the most appropriate kit for your analyser. If more than one tube is suitable, select the smallest volume tube to minimise the blood volume collected.

FBS Kit	Description	Volume $\mu$ l	Tube Length	OD	ID	Heparin type	Concentration IU/ml
R57018-00-SW*	Plastic Capillary, Heparinised	55	59mm +2.0 - 0.0mm	2.05 $\pm$ 0.1mm	1.25 $\pm$ 0.05mm	Na	240
R57024-00-SW*	Plastic Capillary, Heparinised	85	87.0 $\pm$ 0.15mm	2.07 $\pm$ 0.06mm	1.17 $\pm$ 0.02mm	Na	240
R57026-00-SW	Plastic Capillary, Heparinised	100	85.0mm	2.3 $\pm$ 0.05mm	1.20 $\pm$ 0.05	Ca/Na/Li	200

\*Not suitable for lactate measurement on Radiometer analysers as this measurement requires a balanced heparin tube.

### Procedure:

- Using aseptic technique, remove the kit from outer bag, place onto a suitable surface and unfold outer drape to form a sterile field.
- Select a capillary tube from the tube in the kit and load into the holder. Locate the end of the tube into the holder at ❶ and then gently advance into the rear tube grip at ❷. Ensure the tube is secure before sampling. The kit contains a secondary double-ended capillary tube holder if preferred.
- Follow local hospital procedure, prep the perineum and vagina and complete a vaginal examination to establish the cervical dilatation and station of the fetal head.
- Locate the Amnilume™ light module into the slot in the amnioscope and push gently to engage. The light will activate when the light source is correctly fitted. To turn the light off, withdraw the light module.
- Insert amnioscope complete with obturator into the vagina against the fetal presentation ensuring that the cervix is not trapped.
- Remove obturator. Position the bevelled end by rotating endoscope to obtain the best possible seal against the fetal presentation. Rotate to finally position the light source module in the upper quadrant.
- Maintain close contact between amnioscope and scalp to prevent leakage of liquor.
- Clean puncture area with the green x-ray detectable gauze swabs and forceps provided.
- With a clean swab stick, apply a *thin* layer of white soft paraffin (supplied) to puncture area to promote capillary blood aggregation. **Remove excess with a new dry swab stick before skin puncture.**
- The sampling wand comes pre-fitted with an FBS blade. Lightly holding the shaft and rear finger grip (A), push the grip forward to expose the blade (B). **Ensure the blade is locked forward before use.**
- Perform skin puncture through the amnioscope. Either ❶ make a vertical stab, rotate the blade 90° and make a second stab to form a cross or ❷ make a vertical stab, then angle the blade downwards and extend the incision using the front edge of the blade to produce a single incision 4-5mm long.
- Always wait for a good sized (3-4mm) bead to form on the surface before collection. Never slash at the scalp and avoid multiple small stabs.
- To collect the blood sample, depress the button on the wand shaft to retract the FBS blade and angle the capillary tube into the bead.
- Rapidly collect blood sample, filling the tube completely, **without air bubbles**. Collection needs to be expedited to avoid blood clotting. **CAUTION: Do NOT use the sample if the blood clots.**
- Twist lightly to remove the capillary tube from sampling wand and cap both ends with the covers supplied. Take to the analyser **immediately**. Follow the analyser manufacturer's instructions as regards sample mixing. **CAUTION: Avoid distributing air through the sample; samples containing air bubbles may be rejected by the analyser**
- Always ensure the **clean**, non-fetal exposed end is presented to the analyser. Please refer to the analyser instructions for guidance.
- Always ensure that the **correct** measuring mode is selected *before* performing analysis.

**Post sampling:**

18. Haemostasis is accomplished by applying pressure to the puncture site with a swab for one full contraction.
19. If test results so indicate, repeat the procedure. The kit is designed to take up to five individual samples on the same patient.

**Disposal:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment. The Rocket Medical Amnilume™ light module contains lithium-ion based power cells.

This device is not manufactured with natural rubber latex

**Rx** ONLY

STERILE EO

Li-ion

**ROCKET MEDICAL PLC** Sedling Road, Washington,  
England, NE38 9BZ [www.rocketmedical.com](http://www.rocketmedical.com)

EC REP

**Rocket Medical GmbH**  
Am Rosengarten 48,  
15566 Schöneiche.  
Germany

**For Single Patient Use Only.** Do not reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse on another person, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases.

Unless opened or damaged, contents of package are sterile. Store at room temperature. Avoid prolonged exposure to elevated temperatures.

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ZDOCK195 Rev.19 2021-06-17 Copyright© 2010-21 **ROCKET MEDICAL PLC** Todos los derechos reservados. (ES)
1639 ZDOCK195 Rev.19 2021-06-17 Copyright© 2010-21