# THORACIC CATHETERS

## INSTRUCTIONS FOR USE

**Scope:** These instructions cover all Thoracic catheters.

**Device Description:** Thoracic catheter with radiopaque stripe.

**Indications:** For the drainage of air and/or blood and/or pleural fluid and/or pus from the pleural space.

**Contraindications:** Not for non-pleural drainage

**Do not apply more than 10Kpa of thoracic suction to the chest tube.**

**Procedure:**

1. Unwrap catheter, if in an operating theatre use standard protocol for this.
2. Insert the catheter into the pleural space. If in a surgical scenario use a suitable technique. If inserting via blunt dissection use a suitable technique.
3. If using a pull through catheter, insert the catheter into the pleural space by creating a track in the chest wall to pass through. Then pass a suture through the yellow moulded end of the catheter, thread through the track in the chest wall and pull the thoracic catheter into position so the drainage fenestrations remain inside the pleural space, or place forceps through the aforementioned track, grip the yellow moulded end of the catheter and pull into position. Remove the yellow moulded end by cutting it off.
4. Connect to chest drainage collection device.
5. Secure the catheter to chest using a fixation device or suitable suture technique.
6. X-ray patient to ensure correct position of the catheter

**Disposal:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

	This device is not manufactured with natural rubber latex				
					
	<b>ROCKET MEDICAL PLC</b> Sedling Road, Washington, England, NE38 9BZ <a href="http://www.rocketmedical.com">www.rocketmedical.com</a>			<b>Rocket Medical GmbH</b> Am Rosengarten 48, 15566 Schöneiche. Germany	

**CONTINUOUS USE SHOULD NOT EXCEED 28 DAYS.**

**Potential of DEHP Exposure**

Medical procedures using PVC medical devices containing DEHP have the potential to lead to DEHP exposure due to DEHP leaching from the device. The extent of exposure largely depends upon the medical treatments administered, the duration of the treatment, and in the case of plastic blood bags, by the length of storage and the storage temperature. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure. (Full statement at [www.rocketmedical.com](http://www.rocketmedical.com)).

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Rx ONLY

UK  
CA

ROCKET MEDICAL PLC  
Sedling Road, Washington,  
England, NE38 9BZ  
[www.rocketmedical.com](http://www.rocketmedical.com)

Rocket Medical GmbH  
Am Rosengarten 48,  
15566 Schöneiche,  
Germany

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