# Spinal Manometer Sets

## INSTRUCTIONS FOR USE

**Scope:** These instructions cover R55990, R55990-ISO-6 Spinal Manometers and R55995-NRFIT Lumbar Puncture Tray. ISO-6 and NRFIT® versions comply with ISO 80369-6 (as amended).


**Intended Use:** For the measurement of cerebrospinal fluid (CSF) pressure and sampling of CSF to aid the diagnosis of suspected CNS infection, subarachnoid haemorrhage or related neurological condition. This device should only be used by, or under the supervision of, appropriately trained personnel and in conjunction with current national and local clinical practice guidelines. Ref: *NHS National Patient Safety Agency. Safer spinal (intrathecal), epidural and regional devices. Patient Safety Alert, NHS, 2011 and Patient Safety Alert NHS/PSA/RE/2017/004.*

**Contraindications:** This device is **NOT** to be used for the administration of intrathecal agents.

Blind intrathecal puncture is contraindicated in spinal abnormality or injury. Where distortion of normal anatomy is suspected, the use of direct imaging guidance is strongly recommended. Intrathecal puncture is contraindicated in the presence of local skin infections over proposed puncture site, in uncontrolled bleeding conditions or when the patient is receiving anticoagulant therapy.


### Procedure for CSF Pressure Measurement:

1. Position patient either in a curled lateral decubitus position (1) or seated, leaning over a suitable support (2).
2. Locate landmarks: between spinous processes at L4-5, L3-4, or L2-3 levels.
3. Prep and drape the area after identifying landmarks. If required, use suitable local anaesthetic agent to anaesthetise the skin and infiltrate the deeper tissues under the insertion site.
4. Prior to assembly, operate the stopcock 3 or 4 times by rotating it fully 180° to release the seal which may have built up during storage.
5. Assemble the manometer set: locate the lower manometer tube (graduations 3-20cmH<sub>2</sub>O), insert the upper manometer tube (graduations 23-40cmH<sub>2</sub>O) into the luer fitting on the lower tube. Insert the paired manometer tubes into the vertical port of the 3 way tap (as illustrated).

6. To minimise the risk of dural trauma and CSF leakage, the smallest gauge needle that is practical should be used. 25G needles tend to have low flow rates and fill the manometer slowly, therefore 21 or 22G needles are normally recommended for lumbar puncture and CSF collection.
7. Advance the introducer needle (where used) through the skin and fascia, followed by the spinal needle through the deeper tissues to the cord. Typically, a slight pop or reduction in resistance is felt when the dura is punctured.
8. Remove the needle stylet. CSF flow confirms placement into the dural space. Attach the manometer assembly vertically to the needle. Rotate the stopcock to divert the CSF flow from the needle to the manometer whilst closing the proximal port.
9. Allow CSF to flow into the manometer tube, when the fluid level has stabilised, read the value from the manometer column in cmH<sub>2</sub>O. Pressure readings are less reliable if the patient is in the sitting position.
10. Once the CSF pressure reading has been recorded, close the stopcock to the needle and drain the fluid from the manometer tube into collection tubes for subsequent analysis.

R55990 LUER version illustrated


**WARNING: NOT for administration**  
**Do NOT attempt to re-infuse CSF**

11. When sufficient sampling has been completed, remove the spinal needle, manometer assembly and introducer needle (if used) and apply a dressing to the puncture site.
12. Instruct the patient to remain lying down for 1-2 hours before getting up. Lumbar puncture may cause dizziness, profound headache and can disturb balance.

**Disposal:** This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross contamination and infection with transmissible diseases.

This device is not manufactured with natural rubber latex

**R<sub>x</sub>** ONLY

**STERILE**

**ROCKET MEDICAL PLC** Sedling Road, Washington, England, NE38 9BZ [www.rocketmedical.com](http://www.rocketmedical.com)

EC REP

**Rocket Medical GmbH**  
Am Rosengarten 48,  
15566 Schöneiche,  
Germany

**DO NOT RESTERILISE**

Unless opened or damaged, contents of package are sterile.

#### **Kit Contents (R55995-NRFIT)**

- 1 x NRFit® 40cm H2O, stackable spinal manometer (Rocket Medical Plc.)
- 1 x 10ml NRFit® slip syringe (GBUK Healthcare)
- 5 x 10cm x 10cm gauze swabs (Premier Healthcare and Hygiene)
- 1 x 60ml pot (Cups Direct)
- 1 x Skin prep sponge stick (Promepila SA)
- 1 x 2-piece fenestrated drape (Medline Industries Ltd.)
- 1 x square drape 60 x 60 cm (O&M Halyard UK Limited)
- 3 x 30ml Universal Containers with labels (Rocket Medical plc.)
- 1 x 20ml luer lock syringe (BD)
- 1 x Microlance 21G /1.5" green needle (luer) (BD)
- 1 x Surguard 2 25G/38mm orange needle (luer) (Terumo)

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- 1 x 2dílná perforovaná rouška (Medline Industries Ltd.)
- 1 x čtvercová rouška 60 x 60 cm (O&M Halyard UK Limited)