
Rocket® Pudendal Block Needle

INSTRUCTIONS FOR USE

R57040 LUER Fitting Needle illustrated

Scope: These instructions cover: R57040 Rocket® Pudendal Block Needles and derivatives. Sterile, for single use.

Indications: For the introduction of a local anaesthetic agent to create a pudendal nerve block for analgesia for the second stage of labour, repair of an episiotomy or perineal laceration, outlet instrument or for minor surgery of the lower vagina and perineum. This device should only be used by or under the supervision of appropriately trained personnel and in conjunction with current clinical practice guidelines.

Contraindications: Use of the device is contraindicated in the case of the patient's refusal or inability to cooperate; where there is known sensitivity to local anaesthetic agents; in the presence of infection within the ischio-rectal space or the adjacent structures including the vagina or perineum; in patients with coagulation disorders.

Technique:

- Using aseptic technique, remove the needle from the outer bag.
- REMOVE the **transport clip** from the needle.
The depth of submucosal penetration is limited by the device to minimise the risk of injury to the vagina and the fetus.
- The ring clip can be located in either position ① or ② depending on the approach and surgical preference.
- Following local hospital procedure, prep the perineum and vagina using an approved antibacterial solution.
- Palpate the ischial spine either transvaginally or rectally. The ischial spines can be palpated at 4 and 8 o'clock about a finger's length into the vagina. The pudendal nerve runs behind the sacrospinous ligament about 1 cm antero-medial and postero-medial to the ischial spine.
- To perform a left-sided block, palpate the ischial spine with the index finger of the left hand, hold the syringe in the right hand and guide the needle/needle guide between the index and middle finger of the left hand toward the ischial spine

Triple Injection Technique:

- Place the end of the guide beneath the tip of the ischial spine and guide the needle into the vaginal mucosa.
- Aspirate to ensure that the injection is not intravascular.
- Raise a mucosal wheal with local anesthetic.
- Advance the needle through the vaginal mucosa until it touches the sacrospinous ligament, which is located approximately 1 cm medial and posterior to the ischial spine.
- Infiltrate the tissue with local anesthetic.
- Advance the needle further through the sacrospinous ligament for a distance of 1 cm until a loss of resistance is detected. The needle tip now lies in the area of the pudendal nerve.

CAUTION: Care must be taken to avoid intravascular administration into the pudendal vessels which lie lateral to the pudendal nerve.

- Aspirate to confirm the needle placement is not intravascular prior to injection.
- Inject local anesthetic solution into this region.
- Subsequently, withdraw the needle into the guide and move the tip of the guide to just above the ischial spine, reinsert the needle through the mucosa and inject local anesthetic.

Single Injection Technique:

- An alternative technique to the triple injection technique uses a single injection of a larger volume of local anaesthetic. When using this technique omit steps 12-15 above.
- A second alternative single injection technique is performed as above (omitting steps 12-15) but the injection is made 1cm anterior and medial to the ischial spine rather than posterior and medial.
- To perform a right-sided block, repeat these steps holding the syringe in the left hand and guiding the needle/needle guide between the index and middle finger of the right hand toward the ischial spine.
- The patient must be checked bilaterally to establish the extent and efficacy of the block before proceeding with the surgical procedure. The block typically becomes effective about 10 minutes after infiltration.
- Appropriate post procedural monitoring of the patient and the fetus is mandatory. Ensure that provision for resuscitation is made, should adverse reactions or toxicity occur. Exercise caution in the combination and dosage of sedatives and narcotics to avoid possible airway obstruction or aspiration.
- When an early pudendal block is given, the block may later be repeated as long as the maximum dose of local anaesthetic is not exceeded.

(GB)

Disposal

This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

**CAUTION: Do NOT attempt to reattach the transport clip prior to disposal.
Discard the removable finger ring and place the needle and sleeve assembly directly into an approved sharps container.**

This device is not manufactured with natural rubber latex

STERILE EO

Rx ONLY

MD

UK
CA

ROCKET MEDICAL PLC
Sedling Road, Washington,
England, NE38 9BZ
www.rocketmedical.com

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany


For Single Use Only. Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases. **Unless opened or damaged, contents of package are sterile.** Store at room temperature. Avoid prolonged exposure to elevated temperatures.

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Copyright© 2013-2021
All Rights Reserved (GB)
ZDOCK225 Rév. 11 2021-10-20 Copyright© 2013-2021