Hysterocone™ - 6Fg

R57430-CONE

INSTRUCTIONS FOR USE


Intended Use: For the administration and retention of aqueous and oil-based x-ray and ultrasound contrast media within the uterus to demonstrate the uterine cavity and tubal structures during HSG and HyCoSy procedures.

Indications: To aid the diagnosis of fallopian tube blockage, uterine abnormality and other pathological conditions associated with primary and secondary infertility. This device is only for use by or under the supervision of trained personnel and in accordance with local and national guidelines such as those published by the Royal College of Obstetricians & Gynaecologists and the Royal College of Radiologists.

Contraindications: Not for use in patients with a vaginal infection, sexually-transmitted disease or suspected pregnancy.

Preparation:

1. This diagnostic procedure can normally be carried out without medication, although light IV sedation and/or analgesia may be necessary in anxious patients. Use aseptic technique throughout the procedure.
2. Place the patient into the supine position with knees raised, expose the cervix with a Cusco or Sims type speculum. Assess the orientation of the cervix and establish the position of the external cervical os. Swab away any excess mucus with sterile normal saline or an aqueous antiseptic solution as appropriate.
3. Using aseptic technique, open the package and remove the device.
4. Select the desired length of insertion for the intrauterine portion of the cannula by sliding the silicone cone to the appropriate 1cm marking on the cannula to meet the individual patient's needs.
5. Pre-fill the catheter with saline or contrast media to avoid the risk of erroneous diagnosis caused by air bubbles. To aid insertion, a small quantity of sterile water-soluble lubricant may be applied to the tip
6. Locate the distal tip of the Rocket Hysterocone™ cannula into the external cervical os. Using a very gentle rotating motion, advance the cannula into the uterus. The patient may report slight discomfort when the catheter tip passes through the internal os and into the uterine cavity.
7. Locate the silicone cone against the external os to create a fluid tight seal. The inner support tube can be left in place to retain the rigidity of the device or if removed, retain the cone in place with forceps if required.

Warning: when using oil-based contrast media, the cone may become lubricated, sliding more easily on the catheter. To prevent movement during the procedure, place forceps immediately behind the cone. Failure to follow this instruction may result in uterine perforation.

8. Remove the rear cap, attach a 10 or 20ml syringe to the rear port and gently inject sufficient contrast media to adequately demonstrate the uterine cavity and Fallopian tubes.
9. Further volumes of contrast media may be added by turning the tap through 90° to the 'Off' position prior to replenishing the syringe.

To remove the device: Gently withdraw the cannula from the uterus.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health and safety and care of the environment.

This device is not manufactured with natural rubber latex

Rx ONLY

STERILE EO

ROCKET MEDICAL PLC
Sedling Road, Washington,
England, NE38 9BZ
www.rocketmedical.com

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

Unless opened or damaged, contents of package are sterile.

Store at room temperature. Avoid prolonged exposure to elevated temperatures Hysterocone™ is a trademark and Rocket a registered trademark of **ROCKET MEDICAL PLC**

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