Rocket

Pleural Vent Pneumothorax Device

INSTRUCTIONS FOR USE

Scope: These instructions cover all Rocket R54565 Pleural Vent Pneumothorax Device Sets and derivatives. This device should only be used by, or under the supervision of, appropriately trained personnel and in conjunction with current local clinical practice guidelines


Device Description: The Rocket R54565 Pleural Vent Drainage Set contains an aspiration line, insertion pack, and pleural vent assembly consisting of a Veress Needle mounted 8Fg catheter with 6 fenestrations, 1 cm markings commencing 1cm from the distal fenestration, a self-contained one-way valve and vent, an indication diaphragm, a hydrocolloid fixation device, and a needleless port.

The blue indication diaphragm deflects upwards when the catheter tip enters the pleural space and continues to fluctuate with respiration as the intrapleural pressure varies from $> +2\text{mm Hg}$ to $< -2\text{mm Hg}$. This kit does not contain the occlusion cap, R54565-CP, to be ordered separately.

Indications: For the percutaneous introduction of a catheter into the chest for the drainage of air.

Contraindications: Not intended for the drainage of fluid, blood or empyema.



CAUTION: If used for secondary pneumothorax, the patient should be closely monitored due to the possibility of related conditions such as increased fibrin in the pleural fluid, which may lead to the device not performing as intended.


WARNING Not for use with airleaks over 4L/min in volume

Procedure:

1. Ensure that adequate imaging has been performed to confirm the presence of a pneumothorax prior to the use of this drainage device.
2. The recommended position for device use is the 2nd intercostal space in the midclavicular line.
3. Following local hospital policy, prepare the catheter insertion site with an approved solution and drape as required to maintain aseptic technique.
4. Administer appropriate and adequate local anaesthetic to the catheter insertion site and the underlying tissue.
5. Using the scalpel make a skin incision of 4-5mm.
6. Begin to slowly insert the needle through the skin incision over the superior border of the rib and into the pleural space.
7. During insertion through the intercostal muscle, the needle hub window will show RED to indicate the needle tip is exposed.


8. As the needle passes through the parietal pleura a distinct 'click' can be heard as the spring-loaded obturator snaps forward to aid in protecting the internal organs from the needle point.
9. Check the hub window is showing GREEN to indicate the obturator is fully forward.

10. A RED indicator showing or partially showing indicates that the obturator is retracted and the needle tip may be exposed. In this condition, remove the needle, ensure that there is no tissue obstructing the free movement of the obturator and repeat the insertion procedure.


CAUTION: DO NOT ATTEMPT TO REINSERT THE PUNCTURE NEEDLE THROUGH THE DEVICE

As you remove the needle, the blue indication diaphragm will deflect upwards. It then continues to fluctuate with respiration as the intrapleural pressure varies from $> +2\text{mm Hg}$ to $< -2\text{mm Hg}$.

WARNING: Do not over-insert the needle into the chest. Only insert the needle sufficiently to be able to aspirate air.

11. The aspiration of air should be used to verify correct position of the device.
12. When the placement of the needle and catheter has been correctly verified as being in the pleural cavity, remove the needle and syringe from the one-way valve by gently rotating the connection cap anti clockwise. Slowly advance the catheter into the pleura whilst removing the needle from the catheter.
13. When the Veress needle is fully removed, and the catheter is in the pleura, ensure the stitch plate is evenly against the patient's chest wall. At this point, if required, secure the stitch plate to the chest wall using a suture or tape.
14. Attach the valve cap.
15. If fluid drains into the Pleural Vent, this can be removed by attaching a luer lock syringe to the needleless port and withdrawing the fluid.
16. When the indication diaphragm remains in the down position and the patient is asymptomatic the pneumothorax may be resolved. This must be confirmed with imaging before removal of the Pleural Vent.
17. If confirmation of an ongoing airleak is required, the valve system of the Pleural Vent can be occluded using the **R54565-CP** red occlusion cap. This has the same effect of clamping a chest drain. The red occlusion cap should only be used if / as determined by the Clinician reviewing the Patient and should only be left in place under the instructions provided by the Clinician.

WARNING: IF PNEUMOTHORAX HAS NOT RESOLVED AND THE PLEURAL VENT IS TO BE LEFT IN POSITION. **DO NOT LEAVE OCCLUSION CAP IN SITU. REMOVE AS SOON AS POSSIBLE.**
THE PATIENT IS AT RISK OF AN INDUCED TENSION PNEUMOTHORAX.

18. During aspiration of air, if connection to an under-water seal chest drain or suction is required this can be done using R54565-AL.

WARNING: IF ASPIRATION LINE 3 WAY TAP IS LEFT IN THE OFF POSITION, THIS WILL PREVENT THE FLOW OF AIR THROUGH THE UNIT. IF THE PNEUMOTHORAX HAS NOT RESOLVED, THE PATIENT IS AT RISK OF AN INDUCED TENSION PNEUMOTHORAX

DO NOT LEAVE ASPIRATION LINE IN SITU. REMOVE AS SOON AS POSSIBLE.

Reordering Codes	
Rocket Medical Pleural Vent	R54565
Rocket Medical Pleural Vent Occlusion Cap	R54565-CP
Rocket Medical Pleural Vent Aspiration Line	R54565-AL

How to use the Swann Morton Safety Scalpel

1. Grasp the scalpel and carefully extend the blade by moving the slide toward the tip of the scalpel, using the thumb of the hand holding the scalpel.
2. Extend the slider until you reach the positive stop, the slider will fit in the notch when it is completely extended.
3. To retract the blade, grasp the scalpel carefully and move the slide toward the back of the scalpel, using the hand holding the scalpel.
4. You should feel clicks as the blade is retracted and a positive stop once the blade is completely retracted.
5. To retract the blade permanently, move the slider past the notch at the back of the scalpel.
6. Dispose of the Swann Morton Retractable Disposable Scalpel in a puncture resistant container approved for sharps disposal, or in a manner consistent with your hospital procedure.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

	<p>This device is not Manufactured with natural rubber latex</p>									
---	--	---	---	---	---	---	---	---	---	--

	<p>ROCKET MEDICAL PLC Sedling Road, Washington, England, NE38 9BZ www.rocketmedical.com</p>		<p>Rocket Medical GmbH Am Rosengarten 48, 15566 Schöneiche. Germany</p>
---	---	--	--

Potential of DEHP Exposure

Medical procedures using PVC medical devices containing DEHP have the potential to lead to DEHP exposure due to DEHP leaching from the device. The extent of exposure largely depends upon the medical treatments administered, the duration of the treatment, and in the case of plastic blood bags, by the length of storage and the storage temperature. Lack of research in humans means it is difficult to predict the adverse effects of DEHP because certain animal models may not apply to humans. Most studies utilise mice and rats, however breakdown of DEHP in the human body differs to the mechanism observed in these animals. These effects are also seen only at levels far in excess of normal human exposure.

For full statement visit www.rocketmedical.com

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