Rocket® Intrauterine Insemination Catheters:

INSTRUCTIONS FOR USE

R57625 Rocket® DUO™ IUI Catheter illustrated

SCOPE:

These instructions apply to all Rocket Intrauterine Insemination (IUI) Catheters including R57622, R57625, R57625-00-23 and derivatives.

INDICATIONS:

Rocket IUI Catheters: For the placement of human semen into the uterine cavity. This device may only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines and regulations such as those published by the National Institute of Clinical Excellence and the Human Fertilisation & Embryology Authority,

CONTRAINDICATIONS:

Not for use in patients with a vaginal infection, a sexually transmitted disease or suspected pregnancy.

INSTRUCTIONS FOR USE:

Cervical Assessment: The position and patency of the cervix and canal should be assessed 2-3 weeks before the date of initial insemination. If more than 3 months have elapsed since the previous insemination or assessment, then a further examination should be considered. The total uterine length, including endocervix, should be established by uterine sounding or ultrasonography at the initial patient consultation. Patients who have a fully obstructed cervix or impassable endocervical canal should be considered for cervical surgery or alternative fertility treatment before further attempts at IUI are made.

Insemination:

1. Place the patient in the lithotomy position and expose the cervix using a conventional Cusco or Sims speculum. The procedure can normally be carried out without medication, although light sedation may be necessary in very anxious patients.
2. Assess the position of the cervix, establish the position of the external cervical os and swab away any excess mucus with normal saline solution.
3. R57622 is a single part device – the inner catheter cannot be detached. R57625 and R57625-00-23 are two-part devices. Do not remove the inner catheter from the outer sheath as this provides additional support and control.
4. Using aseptic technique, remove the device from the peel pouch and attach a sterile microsyringe.
5. Draw up 1ml of media through catheter into the syringe and discard. Do not flush media from a syringe into the catheter.
6. Draw up nominally 1ml of prepared washed semen into the syringe and catheter assembly.
7. Using aseptic technique, follow one of the two techniques below:
 - a) **No significant anteversion or retroversion, cervical os fully visible:** Locate the distal tip of the IUI catheter into the external cervical os and with a very gentle rotating motion, advance the catheter into the uterus to the required depth. The patient may report slight discomfort when the catheter passes through the internal os and into the endometrial cavity.
 - b) **Anteversion, retroversion with cervical os not fully visible or previous catheter insertion difficult or obstructed:**
 - a. **R57622 Rocket EDL IUI Catheters** contain an internal support wire which aids shaping of the inner catheter
 - b. **R57625 and R57625-00-23 Rocket DUO IUI Catheters** are 2part devices and angulation of the outer sheath can be obtained by gently curving between gloved forefinger and thumb.
8. Once in position, inject the required volume of prepared semen.
9. Whilst maintaining pressure on the syringe plunger, gently withdraw the catheter from the uterus
10. The patient may then be returned to the ward for subsequent discharge.

DISPOSAL: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

ROCKET MEDICAL PLC
Sedling Road, Washington,
England, NE38 9BZ
www.rocketmedical.com

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
ROCKET MEDICAL PLC
Sedling Road, Washington,
England NE38 9BZ
www.rocketmedical.com

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Copyright© 2014-2022. **ROCKET MEDICAL PLC**
Rocket Medical GmbH