

Rocket[®]

R54500-RFE Chest Drain Filter & Coupling Gasket

INSTRUCTIONS FOR USE

Scope: These instructions cover **R54500-RFE**.

Device Description: Individually packed, supplied non-sterile and single use.

Indications: For the single patient use of connecting a Intersurgical breathing filter to R54500/R54567 Family of Chest Drain Bottles.

Contraindications: Not for the connection of a non Intersurgical breathing filter or other apparatus to R54500/R54567 Family of Chest Drain Bottles.

Warning: **BOTH the red coupling and filter require changing every 24 hours.**

Procedure:

1. Place the red coupling gasket on the exhaust port of the R54500/R54567 Family of Chest Drain Bottles.
2. Orientate the filter so the the correct aperture will fit and tightly grip the R54500-RFE.
3. Place the filter on to the coupling gasket, push to secure the fit.
4. If the chest drain bottle is changed, replace the Chest Drain Filter & Coupling Gasket.
5. Replace the filter and red coupling gasket every 24 hours as per Intersurgical instructions.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.



Rx ONLY



ROCKET MEDICAL PLC Sedling Road, Washington,
England, NE38 9BZ www.rocketmedical.com

EC REP

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche,
Germany

NON-STERILE

For Single Use Only. Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases



