HyView[®] HSG and HyCoSy Catheter

INSTRUCTIONS FOR USE

Catheter Product Codes: R57432, R57432-00-ST, R57433, R57433-00-ST.

Intended Use: For the instillation of aqueous or oil-based x-ray and ultrasound contrast media during HSG and HyCoSy procedures.

Indications: To aid in the diagnosis of fallopian tube blockage, uterine abnormality and other pathological conditions associated with primary and secondary infertility. This device is only for use by or under the supervision of trained personnel and in accordance with local and national guidelines.

Contraindications: Not for use in patients with a vaginal infection, sexually transmitted disease, active vaginal bleeding or suspected or confirmed pregnancy.

Duration: Transient, up to 60 minutes.

WARNING: Maximum balloon volume = 1.5 ml.

Do not overfill the balloon as this may cause severe patient discomfort, or balloon rupture and loss of the contrast seal.

Do not use air to inflate the balloon as the balloon material is gas permeable. Use of air to inflate the balloon may lead to a reduction in balloon size over the course of the procedure, leading to loss of the contrast seal.

CAUTION: NOT for use with pressure injectors.

CAUTION: The balloon should be protected from light until use.

Preparation:

1. Use aseptic technique throughout. Follow local hospital policy regarding PPE, which should include provision for the safe handling of contrast media. This diagnostic procedure can normally be carried out without medication, although light IV sedation and/or analgesia may be necessary in very anxious patients.
2. Open the package and remove the device and syringe. Fill the syringe with 1.5ml saline and connect the syringe to the shorter side arm of the catheter (balloon port). Inject 1.5ml of saline into the balloon and lightly palpate with the fingers to confirm balloon patency. Withdraw **all** residual saline from the balloon.
3. Place the patient into the lithotomy position with the knees raised, expose the cervix with a Cusco or Sims type speculum. Assess the orientation of the cervix and establish the position of the external cervical os. Swab away any excess mucus with sterile normal saline or an aqueous antiseptic solution.
4. Pre-form the catheter with the wire stylet in place (applies only to R57432-00-ST and R57433-00-ST versions with integral forming wire).
5. Pre-fill the catheter with contrast media to avoid the risk of erroneous diagnosis caused by air bubbles.
6. To aid insertion, a small quantity of water soluble lubricant may be applied to the catheter tip and balloon.
7. Place the distal tip of the catheter into the external cervical os. Using a very gentle rotating motion, advance the catheter into the uterus, avoiding the uterine wall. The patient may report slight discomfort when the catheter passes through the internal os and into the uterine cavity. The outer sliding sheath may be used to provide extra rigidity to the catheter to aid initial insertion.
8. Using the supplied syringe, fill the balloon with up to **1.5ml of saline**, sufficient to create a seal (Figure. 1). Gently withdraw the catheter until it is located against the internal os. The pressure of the balloon against in the internal os creates a seal and prevents leakage of the contrast media.
9. The catheter may also be inflated in the endocervical canal if insertion beyond the internal os is deemed undesirable.
10. Attach a 10 or 20ml Luer Lock syringe to the rear port, ensuring the tap is in-line with the catheter. Slowly instil 10-20ml of contrast media, sufficient to adequately demonstrate the uterine cavity and fallopian tubes. Further volumes may be added by turning the tap through 90° to the 'Off' position and replenishing the syringe.



Figure 1

To remove the device: Attach the supplied syringe to the balloon port, open the tap and withdraw the plunger to remove all saline from the balloon. Deflation may be slow in devices with stylets. Gently withdraw the catheter from the uterus.

Disposal: This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

ROCKET MEDICAL PLC
Sedling Road, Washington,
England, NE38 9BZ
www.rocketmedical.com

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

DO NOT RESTERILISE

For Single Use Only. Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases. **Unless opened or damaged, contents of package are sterile.**

Distributed in the USA by **ROCKET MEDICAL**

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Copyright© 2022
All Rights Reserved. (en)