# Amniotomy Hook

## INSTRUCTIONS FOR USE

**Scope:** These Instructions apply to the Rocket Medical Amniotomy Hook R57065

**Indications:** Device for the artificial rupture of membranes (ARM) to reduce the duration of labour. This device should only be used by or under the supervision of trained personnel and in conjunction with current clinical practice guidelines.

**Contraindications:** This procedure is contraindicated in placenta previa, vasa previa, in suspected velamentous insertion of the umbilical cord, where the cord is palpable below or near the fetal head, in the presence of malpresentation e.g. transverse lie, in the presence of genital infections (e.g. herpes, Group B streptococcus, gonorrhoea, HIV/AIDS and Hepatitis B), or where the mother is a confirmed carrier of haemophilia and the fetus is either affected or of unknown status. If the fetal head is not fully engaged or the presentation is not determined there is an increased risk of cord prolapse. Ensure that there are no contraindications for vaginal delivery.

### Preparation:

1. Using aseptic technique, remove the hook from the outer bag.
2. Follow local hospital procedure to prepare the perineum and vagina using an approved antibacterial solution.
3. Perform an abdominal palpation prior to verify that the fetus is engaged.
4. Ensure that the mother's bladder is empty.
5. Establish the Fetal Heart Rate (FHR).
6. Position the woman in a modified dorsal position (use a wedge if required), or lithotomy position.
7. Supine hypotension may occur if the woman is asked to lie on her back with the gravid uterus with potential to occlude the inferior vena cava and aorta.
8. Perform a vaginal examination to confirm the fetal presentation and station.
9. Determine the position of the fetus. Exclude any contraindication to performing an ARM.
10. Identify the cervical os and membranes.
11. Insert the index and middle fingers of the dominant hand into the vagina. Using the other hand to hold the device handle, slide the hook between the two internal fingers. Take care to protect the maternal tissues from the hook.
12. When the membranes can be palpated, use the external hand to advance the hook tip to pierce the membranes. Using the external hand, pull the hook back through the membranes to create a tear.
13. After one or two well applied movements of the hook check for amniotic fluid. Do not excessively scratch. Take care to protect the fetus from the action of the hook.
14. Do not immediately withdraw the hand as the risk of cord prolapse is highest during fluid drainage. Observe the amniotic fluid for blood, staining or meconium. Once fluid has drained, verify that cord prolapse has not occurred and withdraw the hand.

**CAUTION: Amniotomy may result in cord prolapse and fetal lacerations.**

### Post artificial rupture of membranes:

1. Observe patient and fetus for signs of distress.
2. Verify that no fetal contusions have occurred during the amniotomy procedure.
3. At intervals, monitor maternal temperature to identify any signs of infection.

### Disposal


This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

**ROCKET MEDICAL PLC**  
Sedling Road, Washington,  
England, NE38 9BZ  
[www.rocketmedical.com](http://www.rocketmedical.com)

**Rocket Medical GmbH**  
Am Rosengarten 48,  
15566 Schöneiche.  
Germany

For Single Use Only. Do not reuse on another person, reprocess or re-sterilise as doing so may compromise the structural integrity of the device, leading to device failure; potentially the cause of serious harm to patients and users. Reuse, reprocessing or re-sterilisation may also result in serious harm to patients and users from cross-contamination and infection with transmissible diseases. **Unless opened or damaged, contents of package are sterile.** Store at room temperature. Avoid prolonged exposure to elevated temperatures.

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