

Rocket

Bone Marrow Aspiration & Biospy Kit

Rocket Bone Marrow Aspiration & Biospy Kits are the direct result of recent user research. The new design includes all those items most commonly identified as required in the ideal ward or department procedure kit.

- **TWO COMPREHENSIVE KITS:** designed to contain the most commonly required items for fast, reliable bone marrow aspiration and marrow biopsy.
- **REDUCES RISK AND STRESS** caused by searching for common items in the ward or department and struggling with missing or worn instruments which are frequently common in reusable trays.
- **REDUCES WASTE:** Combining components into a single kit significantly reduces wasted packaging and reduces wasted time resulting in manufacturing cost savings which means a better utilisation of Trust resources.



- **KIT CONTENTS:**

1 x Rocket T Handle Bone Marrow Aspiration Needle 15G (illustrated)
or
1 x Rocket Marrowgauge Bone Marrow Aspiration Needle 15G
1 x Rocket Trephine Bone Marrow Biospy Needle 11G
3 x 30ml Universal Containers, capped with labels
10 x VFM microscope slides, twinfrost, square edge
2 x 5 slot microscope slide transport boxes
3 x 2ml Syringe luer lock
1 x 12ml Syringe luer lock
1 x 20ml Syringe luer lock

1 x Microlance 21g /1.5" green needle
1 x Microlance 25g 5/8" yellow needle
5 x 10cm x 10cm gauze swabs
1 x 60ml gallipot
1 x 8" skin prep sponge stick
1 x fenestrated drape
1 x square drape 1m²
1 x 5in sharp/blunt scissors
1 x Mepore wound dressing 110mm x 150mm

Ordering Information	Code
Rocket Bone Marrow Biopsy Kits with 'T' Handle Bone Marrow Aspiration Needle Packed in ward procedure tray with sterile drape, in cartons of 5 kits sterile, for single patient use	R56750
Rocket Bone Marrow Biopsy Kits with traditional Marrowgauge Aspiration Needle Packed in ward procedure tray with sterile drape, in cartons of 5 kits sterile, for single patient use	R56751

