


Rocket[®] Marrowgauge[™] Needle



Rocket[®] Marrowgauge[™] Needles provide a highly cost effective solution to the many problems associated reprocessing of conventional marrow aspiration needles.

- **TWO PRIMARY NEEDLE LENGTHS:** to meet a wide range of patient needs.
- **CONTROLLED DEPTH** of penetration with 2 sleeves allowing sternal and iliac crest puncture with total control of penetration depth.
- **GROUND TAPER TIP** with precision ground tip to enable easy penetration of bone.
- **ERGONOMIC HANDLE** sits comfortably in the operator's hand during initial penetration.
- **SIMPLE** locking mechanism secures stylet to ensure clean controlled bone penetration.

Ordering Information	Sternal Iliac Crest	NHS Order Code	Code
Rocket[®] Marrowgauge[™] Needles 15g, with moulded body, retaining cap, stylet and removable sleeves. Peel pouch packed in cartons of 10 needle sets, sterile, for single patient use	10-22mm 22-39mm	FTR324	R56780
	35-47mm 47-64mm		R56781

