

# Rocket NPSA

## Spinal Manometers & LP Procedure Packs


Produced in partnership with the manufacturers of BD UniVia® and SURETY® Neuraxial Connectors, Rocket Medical have developed a range of Spinal Manometer sets and Lumbar Puncture Procedure Pack specifically designed to meet the improved safety requirements of the NPSA Alert: *Safer Spinal (intrathecal) epidural and regional devices. Parts A & B. (as amended).*

- **REDUCES RISK:** Combining all the necessary compatible components into a single kit significantly reduces the risk of mismatch or accidental connection to non-approved systems.
- **CLEAR, STURDY CLEAR CONSTRUCTION** for reliable, secure assembly and ease of use. Ready for immediate use for when time is precious.
- **TWO PART MAMOMETER ASSEMBLY** with extension allows CSF pressure monitoring between 3-55cmH<sub>2</sub>O. Supplied with 3 way stopcock for attachment to a BD UniVia® or SURETY® compatible NPSA compliant spinal needle and syringe.
- **ALSO AVAILABLE IN A COMPREHENSIVE KIT:** designed to contain the most commonly required items for fast, reliable lumbar puncture, CSF Pressure Measurement and CSF Sampling in the department or ward setting. (SURETY® only)
- **PACK CONTENTS:**



- 1 x SURETY® compatible 30cm Spinal Manometer Set
- 1 x SURETY® compatible 22G x 9cm Quinke spinal needle
- 1 x 10ml SURETY® compatible syringe
- 5 x 10cm x 10cm gauze swabs
- 1 x 60ml gallipot
- 1 x 8" skin prep sponge stick

- 1 x fenestrated drape
- 1 x square drape 1m<sup>2</sup>
- 3 x 30ml Universal Containers with labels
- 1 x 20ml Syringe luer lock
- 1 x Microlance 21G /1.5" green needle (luer)
- 1 x Microlance 25G/5/8" yellow needle (luer)



Ordering Information	Code	NHS SC Code
<b>Rocket BD UniVia® NPSA Spinal Manometer Set</b> 1 x Rocket Spinal Manometer Set with 3 way stopcock and extension for single patient use, packed in cartons 10 units	<b>R55990-BD</b>	
<b>Rocket SURETY® NPSA Spinal Manometer Set</b> 1 x Rocket Spinal Manometer Set with 3 way stopcock and extension for single patient use, packed in cartons 10 units	<b>R55990-VN</b>	<b>FTP135</b>
<b>Rocket SURETY® NPSA Lumbar Puncture Procedure Kit</b> complete ward or department procedure kit for single patient use, packed in cartons 5 units	<b>R55995-VN</b>	

Rocket and the Rocketmedical logo are registered trademarks of Rocket Medical plc. Surety® is a registered trademark of Intervene Ltd. UniVia® is a registered trademark of BD

