

Rocket®

Bone Trephine™ Needle



The Rocket® Bone Trephine™ Needle is designed to give safe, reliable penetration performance combined with excellent comfort and sample retention characteristics.

Each set features:

- ◆ **INNER STYLET** with precision ground tip to enable easy penetration of bone.
- ◆ **ERGONOMIC HANDLE** sits comfortably in the operators hand during initial penetration.
- ◆ **COVER PLATE** to allow controlled collection of sample after stylet has been removed.
- ◆ **SIMPLE** locking mechanism secures stylet to ensure clean controlled bone penetration.
- ◆ **GROUND TAPER TIP** on outer needle ensures reliable specimen retention.
- ◆ **SAMPLE EJECTOR** supplied to quickly release specimen without damage and to minimise the risk of needle stick injury
- ◆ **SAFETY COVER** to protect the potentially contaminated needle point during sample removal for enhanced user safety.

Ordering Information	Size	NHS Logistics	Code
Rocket® Bone Trephine Needle Set: Needle, trocar, blank cover plate, sample ejector, safety cover. Sterile, peel pouch packed in cartons of 10 needle sets, for single patient use Sizes: 7g: 4.45mm OD x 100mm, 11g: 2.95mm OD x 100mm, 14g: 2.05mm OD x 40mm - Paediatric	7FG		R56791
	11FG	FTR349	R56790
	14FG Paediatric		R56792

