


Rocket Rolonspec™

Disposable Vaginal Speculum

The Rocket RolonSpec™ combines cost effective patient comfort with excellent access and vision

- **HIGH QUALITY MOULDING** combines blade strength with smooth edge profiles for maximum access and patient comfort.
- **SILENT RATCHET** operation. Single handed operation allows the speculum to be opened without engaging the ratchet until full distension has been achieved – reducing patient stress and discomfort.
- **INDIVIDUALLY WRAPPED**, sterile or non-sterile in cartons of 25 units.



Ordering Information	Code
Rocket Rolonspec™ Disposable Vaginal Cusco Speculum in cartons of 25 units, sterile, for single patient use	R50600
Rocket Rolonspec™ Disposable Vaginal Cusco Speculum in cartons of 25 units, non-sterile, for single patient use	R50601

