

# CUSTOMER INFORMATION BULLETIN


## Rocket HSG, HyCoSy & Sonohystogram Catheters

Dear Customer,

With immediate effect, all Rocket HSG Catheters are fully available directly from our preferred logistics partner CH2. These product were formally distributed by Gytech Pty.

Rocket HSG, HyCoSy & Sonohystogram Catheters still have these well-established and trusted features:

- **LATEX FREE** balloon minimises risk of patient tissue reactions.
- **REMOVABLE STIFFENING SHEATH** to aid insertion in the heavily retroverted or antiverted cervix.
- **CONTRAST PATH TAP** ensures retention of contrast when further volumes of media are required.
- **SMOOTH TIP & LOW PROFILE BALLOON** ensures low trauma insertion and removal – ideal for the nulliparous cervix.
- **INTEGRAL WIRE STYLET** for improved malleability without using tiresome removable stylets
  - Improves ease of insertion
  - Increased value for money
- **COMPATIBLE** with all aqueous based contrast media.
- **AVAILABLE FROM OUR PREFERRED AUSTRALIA LOGISITICS PARTNER**



CLIFFORDHALLAMHEALTHCARE

1300 720 274 www.ch2.net.au

Ordering Information	Size	Code	CH2 Code
<b>Rocket HSG &amp; Sonohystogram Catheter</b> Length 280mm, balloon inflation syringe pre-attached. Maximum balloon volume: 3ml. Sterile, for single patient use, peel pouch packed in cartons of 10 catheters.	5FG	<b>R57422</b>	<b>2067145</b>
	7FG	<b>R57424</b>	<b>2067153</b>
<b>Rocket HSG &amp; Sonohystogram Catheter with integral wire stylet.</b> Length 280mm, balloon inflation syringe pre-attached. Maximum balloon volume: 3ml. Sterile, for single patient use, peel pouch packed in cartons of 10 catheters.	5FG	<b>R57422-00-ST</b>	<b>2067161</b>
	7FG	<b>R57424-00-ST</b>	<b>2067170</b>

Customer Information Bulletins on all of our products can be found on our web site at:  
[http://www.rocketmedical.com/Rocket\\_support\\_CIB\\_main.htm](http://www.rocketmedical.com/Rocket_support_CIB_main.htm)

**If you have any questions, please contact the Customer Services Team for assistance:**  
 Tel: +61 431 146448 (NSW, ACT) +61 404 532746 (Other territories) [aus@rocketmedical.com](mailto:aus@rocketmedical.com)

