

Rocket Filter Sets

Rocket Filter Sets offer a simple, cost effective solution to the contamination risks posed by reusable water traps or multiple use consumable filter systems.

- **REMOVES CONTAMINATION RISK** from irregular or uncontrolled decontamination of reusable bottle sets or multiple use consumable filters.
- **PATIENT TUBE SET:** 2.0m patient tube set with integral luer lock connector to ensure secure, reliable connection to all standard luer-lock fitting aspiration needle sets.
- **HYDROPHOBIC BARRIER** provides pump protection by preventing accidental aspiration of fluid into the device.
- **SUITABLE FOR USE WITH:**
 - Craft™ Suction and Duo Vac
 - *genX Pioneer Pro Pumps.*
 - COOK K-MAR-5100, K-MAR-5200 and derivatives



FILTER SPECIFICATION
 68mm x 64mm
 Effective filtration area = 24.6mm²
 0.027micron,
 BFE 99.9999%*
 VFE 99.9999%*
 High flow /low resistance Insufflation/Gas/Air Filter
 Hydrophobic barrier.
 Inner core: silicon treated paper

Description	Code
Rocket Filter Set: suitable for Craft™ Suction Pumps, <i>genx</i> Pioneer Pro and derivatives, sterile, for single patient use, cartons 10 units	R57678
Rocket Filter Set: suitable for COOK K-MAR-5100, K-MAR-5200 and derivatives, sterile, for single patient use in cartons of 10 units.	R57679

