

# Rocket Uterine Sounds



- **WELL TOLERATED:** Only 3.5mm OD distal bead and 2.5mm distal shaft for easier, low trauma insertion.
- **MALEABLE DISTAL SHAFT:** allows for increased angulation for easier insertion into the retroverted or antverted uterus.
- **GRADUATED:** 1cm depth markers commence 40mm from the distal tip. Proximal marker ring positioned 90mm from distal tip.
- **NEW HANDLE DESIGN:** increases grip and tip control whilst reducing bulk to improve visualisation of the cervix.
- **STERILE:** Supplied sterile, individually wrapped for single patient in cartons 100 units.

Ordering Information	Code
<b>Rocket Uterine sound.</b> Sterile, for single patient use packed in cartons 100 units	<b>R57420</b>

