

# CUSTOMER INFORMATION BULLETIN

## Rocket FBS Amnioscope with Amnilume: R57016-LUME

### Dear Customer,

We are making improvements to our **R57016-LUME Rocket FBS Amnioscope and Amnilume** by introducing a new design which has these benefits:

- **Light module is easier to insert and remove from the amnioscope.**
- **Light module occupies less space inside the scope - giving more room for instrumentation and improved visibility.**
- **More rounded shape to the obturator - improves comfort and ease of insertion.**
- **Retains the same high output LED of the current Rocket Amnilume.**

The changes, which affect both amnioscope and light source modules, are being introduced from **June 2015** onwards and will initially affect the R57016-LUME and the R57016-PLUS amnioscope sets.

The current design (pictured lower right) will be phased out as stocks of existing components are used up. The new design will be progressively introduced across the whole Rocket FBS Kit range over the next 2-3 months.

**Customer Information Bulletins** on all of our products can be found on our web site at:

[www.rocketmedical.com](http://www.rocketmedical.com)

**If you have any questions, please contact the Customer Services Team for assistance:**

**Tel: +44 191 419 6988.**

**Email: [customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com)**



Original amnioscope & light source design, now being replaced