

CUSTOMER INFORMATION BULLETIN

Rocket FBS Kits

Design changes to the Rocket Amnioscope and Amnilume™

Dear Customer,

We are making improvements to the Rocket **Amnioscope and Amnilume™** contained in our **Rocket FBS Kit** and will be introducing a new design which has these benefits:

- **Light module is easier to insert and remove from the amnioscope.**
- **Light module occupies less space inside the scope - giving more room for instrumentation and improved visibility.**
- **More rounded shape to the obturator - improves comfort and ease of insertion.**
- **Retains the same high output LED of the current Rocket Amnilume.**

The changes, which affect both amnioscope and light source modules, are being introduced from **June 8th 2015** onwards and will affect:

All standard Rocket FBS Kit codes: R57012, 013, 018, 024, 026, 027, 028

and all **Rocket FBS Sampling Wand FBS Kits** with codes ending **R570..-00-SW**

The current design (pictured lower right) will be phased out as stocks of existing components are used up.

Unfortunately, original style Amnilume modules are not compatible with the new design Amnioscope. If you have stock of either R57004 Amnilume™ Light modules and/or R57016-PLUS Amnioscopes, please contact Customer Services for advice.

Customer Information Bulletins on all of our products can be found on our web site at: www.rocketmedical.com



Original amnioscope & light source design, now being replaced

If you have any questions, please contact the Customer Services Team for assistance:

Tel: +44 191 419 6988. Email: customerservices@rocketmedical.com