

CUSTOMER INFORMATION BULLETIN

Rocket® R54500 & R54509 Single Chest Drainage Bottle

Dear Customer,

As part of our programme of continuous improvement we are pleased to introduce an update for the R54500 & R54509 Chest Drainage bottles.

- There is now an automatic positive pressure relief valve.
- If the flow of escaping air is blocked, the positive pressure relief valve will allow the air to escape and protect the patient from potential harm.



Positive Pressure Blow Off Valve

Customer Information Bulletins on all of our products can be found on our web site at: www.rocketmedical.com

If you have any questions, please contact the Customer Services Team for assistance:

Tel: +44 191 419 6988. Fax: +44 419 6889. Email: customerservices@rocketmedical.com

