

# CUSTOMER INFORMATION BULLETIN

## Rocket® R57065 Amniotomy Hook - Upgrade



### Dear Customer,

We are making improvements to our **R57065 Rocket Amniotomy Hook** by introducing an entirely new moulding.


The existing manufacturing tools have reached the end of their service life and are being replaced.

As part of our continuous improvement programme, we have taken the opportunity to re-design the overall handle shape and the increase the rigidity of the hook profile.

The illustrations (right) show the change in hook shape.

The changes, which affect only the R57065 Rocket Amniotomy Hook, are being introduced into all markets from **September 1st 2015** onwards.

New tip profile



Illustrations not to scale

Original tip profile



**Customer Information Bulletins** on all of our products can be found on our web site at:

[www.rocketmedical.com](http://www.rocketmedical.com)

If you have any questions, please contact the **Customer Services Team** for assistance:

Tel: +44 191 419 6988.

Email: [customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com)