

CUSTOMER INFORMATION BULLETIN

Dear Customer,

As part of our continuous improvement programme, we are upgrading our **Rocket® FBS Kits** by introducing a new longer, 'double ended' capillary tube holder. Standard Rocket FBS Kits contain one capillary tube holder and Rocket Sampling Wand FBS Kits contain a capillary tube holder in addition to the sampling wand.

This upgrade is in direct response to end user requests for a double ended tube holder to allow users to perform secondary tube fillings without the need to remove and replace the capillary tube as required with the present single ended design.

The increasing usage of our longer 220mm PLUS size amnioscopes has also resulted in requests for a longer tube holder and this we have included in this upgrade by increasing the tube holder length by a 5cm.

These changes are an upgrade to current product and will be phased into future production as existing stocks of kits are delivered. Existing stocks with single ended tubes holders are still fully usable and can continue to be shipped to customers.

The kit compatibility table below has been recently updated to include the Radiometer ABL™ 90 FLEX PLUS.

Compatible Analyser	Recommended sample volume	Rocket® Standard FBS Kit	Rocket® FBS Kit with Sampling Wand
ABL™ 5, 500, 70, 77	65-85µl	R57024	R57024-00-SW
ABL™ 700 Series with tube volume selected	>35µl	R57018	R57018-00-SW
ABL™ 90 FLEX	65µl	R57024	R57024-00-SW
ABL™ 90 FLEX PLUS	45µl	R57018	R57018-00-SW
ABL™ 800 FLEX	>85µl	R57012 R57013 R57024	R57012-00-SW R57013-00-SW R57024-00-SW
ABL™ 800 FLEX with FLEXMODE selected	55µl	R57018	R57018-00-SW
Instrumentation Laboratory (IL) GEM™ 3000 Series	170µl	R57027	N/A
Instrumentation Laboratory (IL) GEM™ 4000 Series	65µl	R57024	R57024-00-SW
ROCHE OMNI c, OMNI s cobas b 121, b123, b221 systems	35µl 85µl	R57018 R57024	R57018-00-SW R57024-00-SW
SIEMENS® (Bayer) Rapidlab™ Series	100µl	R57026	R57026-00-SW
NOVA Biomedical StatStrip® Lactate	6µl	R57028	R57028-00-SW

Customer Information Bulletins on all of our products can be found on our web site at: www.rocketmedical.com

If you have any questions, please contact the Customer Services Team for assistance:
Tel: +44 191 419 6988. Email: customerservices@rocketmedical.com

