

## CUSTOMER INFORMATION BULLETIN

### Dear Customer,

We are making improvements to all of our **Rocket®** and **NuSurgix® FBS Sampling Wands** by introducing a new blade and blade mounting.

The existing manufacturing tools have reached the end of their service life and are being replaced and upgraded.

As part of our continuous improvement programme, we have modified the blade width and added a dual bevel to increase sharpness and improve penetration of the fetal scalp.

The image (right) shows the blade holder of a NuSurgix® Sampling Wand, the new design (top) with the older, outgoing type (bottom).

We retain the 2mm penetration depth of our standard Rocket FBS blade for added safety.

These design changes will also apply to Rocket® FBS Sampling Wands which will now include the new blade in a modified housing.

Functionally the sampling wand remains unchanged.

The changes, which affect all NuSurgix® and Rocket® FBS Sampling Wands and FBS Kits with Sampling Wands will come into effect from **1st November 2015** onwards.

The new design components will be introduced as existing stocks of finished wands and kits are used up.



**IMPORTANT: These changes only affect combined Sampling Wands, they do NOT affect R57010 and R57011 Rocket® Plug-fit and Screw-in type FBS Blades**

**Customer Information Bulletins** on all of our products can be found on our web site at: [www.rocketmedical.com](http://www.rocketmedical.com)

**If you have any questions, please contact the Customer Services Team for assistance:**

**Tel: +44 191 419 6988. Email: [customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com)**