

## Rocket<sup>®</sup> IPC<sup>™</sup> Insertion Kits R54400-I6-40 & R54400-I6-MI

**Dear Customer,**

As part of our programme of continuous improvement we are pleased to introduce an update for the Rocket Medical range of IPC<sup>™</sup> Insertion Kits.

The stiff plastic 'tunneller' has been replaced with a malleable metal tunneller in the product codes above.



**Customer Information Bulletins**  
for all of our products can be found on our web site at:  
[www.rocketmedical.com](http://www.rocketmedical.com)

**If you have any questions, please contact the  
Customer Services Team for assistance:**

**Tel: +44 191 419 6988. Email: [customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com)**