


Locking hub on all R58900 series Rocket® Introducer Needles

Description		Length	Code
<b>Rocket® Introducer Needles</b> 5FG outer sheath, 19G inner cannula tube, 21G tri-facet trocar needle Sterile, per 10	R/O	19.5cm	<b>R58900</b>
	Clear	19.5cm	<b>R58900-00-KL</b>
	R/O	15cm	<b>R58900-00-15</b>
	Clear	15cm	<b>R58900-KL-15</b>
<b>Rocket® DUO™ Introducer Needle</b> 5FG outer sheath, 19G inner chiba needle Sterile, per 10	Clear	13cm	<b>R58900-CH-15</b>
<b>Rocket® Fixed Core Guidewire, PTFE</b> coated, 3mm J tip, dispenser loaded. Sterile, per 10	0.035"	80cm	<b>R58905-35-80</b>
<b>Rocket® Radiopaque Dilators</b> luer lock, polyethylene, barium loaded tapered dilators Sterile, per 10	5-16FG	21cm	<b>R58921-05-21</b> <b>R58921-06-21</b> <b>R58921-07-21</b> <b>R58921-08-21</b> <b>R58921-12-21</b> <b>R58921-13-21</b> <b>R58921-14-21</b>
<b>Rocket® Access Set</b> PTFE J tip wire guidewires: 0.018", 0.035" 4,5,6,7,8FG radiolucent dilators, 18G Introducer needle Sterile, per 10		80cm 21cm 7cm	<b>R58912</b>
<b>Rocket® Drainage Connection Set</b> Funnel end connector, tubing and luer stop- cock with rotating locking ring. Sterile, per 10		30cm	<b>R58915</b>



All Rocket® R58900 series Rocket® Introducer Needles feature locking hubs and silicone lubrication for easier inner needle withdrawal