

# Rocket® Amniocentesis Exchange Set



- ◆ **HIGH RESOLUTION ECHO BEVEL DESIGN** to give smoother, easier tissue penetration with pinpoint ultrasound definition of the needle tip in either amniocentesis and CVS.
- ◆ **NEEDLE MARKINGS** and depth marker help ensure correct positioning and minimise the risk of over insertion.
- ◆ **INTEGRAL TUBING SET** – purpose designed 10cm extension tube set is included in every pack for simple, safe attachment to the needle without the need to find costly off-the-shelf items which cannot offer guaranteed compatibility and security.
- ◆ **10ml 3 RING SYRINGE & 3 WAY STOPCOCK INCLUDED:** - each set contains a 10ml, 3 ring syringe and widebore 3 way stopcock for comfort and security when repeated aspiration/infusion is required.
- ◆ **SECURE LUER LOCK FITTINGS** ensure the needle, extension set and syringes remain correctly and securely connected throughout
- ◆ **CE MARKED, PURPOSE DESIGNED DEVICE**, designed specifically for amniocentesis and transabdominal chorionic villus sampling which avoids the risks of inappropriate use of non specific generic devices.



The Rocket® Amniocentesis Exchange Set has been developed with the assistance of a leading UK Fetal Medicine Unit and is the result of a program of continuous product development and design improvement.

**Description**
**Rocket® Amniocentesis Exchange Set**

echogenic amniocentesis needle with 1cm markings, 10cm luer lock extension tube set, 3 way- stopcock, 3 ring syringe.  
Sterile for single use, packed in cartons of 10 sets

**Order Code**
**R57035**