



# Rocket® Amniotomy Hooks



- ◆ **NEW RE-DESIGNED** shaft and distal hook to improve rigidity and hook profile
- ◆ **GUARDED HOOK:** located at a controlled depth to permit perforation of the membranes without trauma to the surrounding tissue
- ◆ **CURVED:** anatomically shaped to ensure a positive grip in the users hand and to follow the normal curved pelvic outlet. for precise finger location and positioning.
- ◆ **STRAIGHT: SURGI-HOOK™** with contoured handle to indicate angle of perforating point, Hook faces away from direction of insertion to protect against accidental tear.



Rocket Medical plc have a long and enviable reputation for reliability and innovation in the Labour Ward and our acquisition of the Surgicraft® range complements our existing fetal and labour ward products.

As a result, we can now offer two amniotomy hook versions to meet your individual needs

Description	Order Code
<b>Rocket® Amniotomy Hook</b> , curved, 25cm. Sterile, for single use, individually packed, 100 per carton.	<b>R57065</b>
<b>Rocket® Surgi-Hook</b> , straight 26cm. Sterile, for single patient use, individually packed,	<b>R57066</b>