


The Rocket® Band-It™ single use multiband RBL haemorrhoid ligator comes pre-loaded with 5 latex or hypoallergenic ligature bands for single handed deployment of up to 5 bands without the inconvenience of re-loading.

- **READY TO USE:** self-contained, pre-loaded with 5 latex or hypoallergenic haemorrhoid bands
- **COMPLETE SET** supplied with soft, highly flexible suction tube for reduced drag.
- **INDEXED DEPLOYMENT WHEEL** makes band delivery quick and easy to control
- **SIMPLE OPERATION:** Designed for single handed deployment with either left or right hand.

Description	Style	Order Code
Rocket® Band-It™ RBL Haemorrhoid Ligator. Non-sterile, 20 units per carton.	Latex bands with suction tube	R53702 CH2: 2211063
	Hypoallergenic bands with suction tube	R53703 CH2: 2211071


CLIFFORDHALLAMHEALTHCARE

 1300 720 274 www.ch2.net.au