



- ◆ Used in over 90% of UK hospitals.
- ◆ Safe, quick and easy to apply without the need for applicators.
- ◆ Ribbed finger grip for easier placement.
- ◆ Can be applied with as little as 1.5cm cervical dilation for earlier monitoring.
- ◆ Three point contact for excellent trace quality
- ◆ Automatic control of depth of penetration of needle for low scalp trauma.
- ◆ Cannot detach accidentally when correctly applied.
- ◆ Now also available with push connectors to fit the latest fetal monitoring systems.
- ◆ Easy to reposition and remove.



The Rocket Copeland[®] name and Fetal Scalp Monitoring are synonymous and for over 20 years Copeland[®] Fetal Scalp Electrodes have provided safe and reliable recording of the fetal heart rate during labour.

It is the preferred choice when internal fetal scalp monitoring is indicated.

“The Copeland[®] performed the best overall, and particularly in respect to attachment to the scalp.”

On the basis of the results of this trial we feel that the Copeland[®] is the electrode of choice.”

Ref: British Journal of Obstetrics and Gynaecology, April 1992. Vol99 pp302-306

Description	NHS SC	Order Code
Rocket[®] Copeland[®] Fetal Scalp Electrode twin wire type Sterile, for single use, individually packed, 25 per carton.	FDK221	R57008-00-FE
Rocket[®] Copeland[®] Fetal Scalp Electrode, Quik connector Sterile, for single use, individually packed, 25 per carton.	FDK097	R57008-00-CN
Rocket[®] Copeland[®] Fetal Scalp Electrode for Philips Avalon FM30 & FM50 systems Sterile, for single use, individually packed, 25 per carton.	FDK557	R57008-00-FM