



R57026 with sampling wand illustrated

- ◆ **RELIABLE, COMPATIBLE ANALYSIS:** Rocket FBS Kits are specifically designed for use with Radiometer® ABL™, Siemens® Rapidlab, Instrumentation Laboratory® GEM™, ROCHE OMNI, cobas, Nova Biomedical lactate and other analysers
- ◆ **MAXIMUM ILLUMINATION NOW COMES AS STANDARD:**
To give you brighter, more reliable visualisation every kit contains the latest **Rocket Amnilume™ Light Source.**
- ◆ **EASIER SAMPLE COLLECTION:** The Rocket Sampling Wand combines blade and capillary tube holders to help solve the problem of relocating the puncture site when transferring between blade and tube holder. It is as easy to use a ball point pen.
- ◆ **CAPILLARY TUBE DISPENSERS:** R57028-00-SW series kits contains 3 individual capillary tube dispensers. Simply remove the capillary tube, complete with blood sample, from the Sampling Wand or capillary tube holder, insert into the adaptor on the tube dispenser bulb and lightly depress the bulb to deliver the blood sample directly to the test strip



At Rocket Medical we know our customers are our greatest source of new ideas and product improvements.

Our latest FBS Kits incorporate your user feedback and will mean better visualisation, easier scalp puncture and sample collection all combined in a smaller tray to save on packaging and storage space.

CONTENTS:

- Rocket Amnioscope and Amnilume™ Single Use Light Source,
 - Rocket Sampling Wand - combined blade and capillary tube holder or
 - Separate tube holder and blade holder (R57012 only)
- 5 x 10cm x-ray detectable gauze swabs - green.
- 6 swab sticks
- 2.5gms petroleum jelly
- sponge holding forceps
- high heparin capillary tubes
- spare double-ended capillary tube holder
- tube covers
- 3 x capillary tube dispensers (R57028-00-SW only)
- Packed in cartons of 10 units, for single use

Compatible Analyser	Recommended sample volume	Rocket® FBS Kit
ABL™ 800 FLEX	>55µl	R57012
ABL™ 90 FLEX PLUS ABL™ 800 FLEX with FLEXMODE selected ROCHE® OMNI c, OMNI s cobas b 121, 123, 221 systems	35µl - 45µl	R57028-00-SW
SIEMENS® (Bayer) Rapidlab™ Series	100µl	R57026-00-SW
NOVA Biomedical StatStrip Lactate Pro	6µl	R57028-00-SW