



At Rocket Medical we know our customers are our greatest source of new ideas and product improvements.

Our latest FBS Kits incorporate your user feedback and will mean better visualisation, easier scalp puncture and sample collection all combined in a smaller tray to save on packaging and storage space.

CONTENTS:

- Rocket Amnioscope and Amnilume™ Single Use Light Source,
- Rocket Sampling Wand - combined blade and capillary tube holder
- 5 x 10cm x-ray detectable gauze swabs - green.
- 6 swab sticks
- 2.5gms petroleum jelly
- sponge holding forceps
- high heparin capillary tubes
- spare double-ended capillary tube holder
- tube covers
- Packed in cartons of 10 units, for single use

- ◆ **RELIABLE, COMPATIBLE ANALYSIS:** Rocket FBS Kits are specifically designed for use with Radiometer® ABL™, Siemens® Rapidlab, Instrumentation Laboratory® GEM™, ROCHE OMNI, cobas, Nova Biomedical lactate and other analysers
- ◆ **MAXIMUM ILLUMINATION NOW COMES AS STANDARD:**
To give you brighter, more reliable visualisation every kit contains the latest **Rocket Amnilume™ Light Source.**
- ◆ **EASIER SAMPLE COLLECTION:** The Rocket Sampling Wand combines blade and capillary tube holders to help solve the problem of relocating the puncture site when transferring between blade and tube holder. It is as easy to use a ball point pen.

Compatible Analyser	Recommended sample volume	Rocket® FBS Kit with Sampling Wand	NHS Order codes
ABL™ 5, 500, 70, 77, 80 FLEX	65-85µl	R57024-00-SW	FFY056
ABL™ 700 Series with tube volume selected	35µl	R57018-00-SW	FFY055
ABL™ 800 FLEX	>85µl	R57012-00-SW R57013-00-SW R57024-00-SW	FFY053 FFY054 FFY056
ABL™ 90 FLEX	65µl	R57024-00-SW	FFY056
ABL™ 90 FLEX PLUS	45µl	R57018-00-SW	FFY055
ABL™ 800 FLEX with FLEXMODE selected	35µl	R57018-00-SW	FFY055
Werfen - Instrumentation Laboratory GEM™ 4000, 5000 Series	65µl	R57024-00-SW	FFY056
ROCHE® OMNI c, OMNI s cobas b 121, 123, 221 systems	35µl 85µl	R57018-00-SW R57024-00-SW	FFY055 FFY056
SIEMENS® (Bayer) Rapidlab™ Series	100µl	R57026-00-SW	FFY057
NOVA Biomedical StatStrip Lactate	6µl	R57028-00-SW	N/A