



At Rocket Medical we know our customers are our greatest source of new ideas and product improvements.

Our latest FBS Kits incorporate your user feedback and will mean better visualisation, easier scalp puncture and sample collection all combined in a smaller tray to save on packaging and storage space.

- ◆ **RELIABLE, COMPATIBLE ANALYSIS:** Rocket FBS Kits are specifically designed for use with Radiometer® ABL™, Siemens® Rapidlab, Instrumentation Laboratory® GEM™, ROCHE OMNI, cobas, Nova Biomedical lactate and other analysers
- ◆ **MAXIMUM ILLUMINATION NOW COMES AS STANDARD:** To give you brighter, more reliable visualisation every kit contains the latest **Rocket Amnilume™ Light Source.**
- ◆ **COMPREHENSIVE:** each specifically developed entry-level kit has been designed in conjunction with the world's leading manufacturers of blood gas analysers to offer reliable, safe and cost effective fetal blood sampling.

CONTENTS:

- Rocket Amnioscope and Amnilume™ Single Use Light Source,
- Rocket blade holder
- 2 x FBS blades
- double-ended capillary tube holder
- 5 x 10cm x-ray detectable gauze swabs - green.
- 6 swab sticks
- 2.5gms petroleum jelly
- sponge holding forceps
- high heparin capillary tubes
- tube covers
- Packed in cartons 10 units, for single use

Compatible Analyser	Recommended sample volume	Rocket® Standard FBS Kit	NHS Order code
ABL™ 5, 500, 70, 77, 80 FLEX	65-85µl	R57024	FFY079
ABL™ 700 Series with tube volume selected	>35µl	R57018	FFY078
ABL™ 90 FLEX	65µl	R57024	FFY079
ABL™ 90 FLEX PLUS	45µl	R57018	FFY078
ABL™ 800 FLEX	>85µl	R57012 R57013 R57024	FFY076 FFY077 FFY079
ABL™ 800 FLEX with FLEXMODE selected	55µl	R57018	FFY078
Werfen - Instrumentation Laboratory GEM™ 3000 Series	170µl	R57027	N/A
Werfen - Instrumentation Laboratory GEM™ 4000, 5000 Series	65µl	R57024	FFY079
ROCHE OMNI c, OMNI s cobas b 121, b123, b221 systems	35µl 85µl	R57018 R57024	FFY078 FFY079
SIEMENS® (Bayer) Rapidlab™ Series	100µl	R57026	N/A
NOVA Biomedical StatStrip® Lactate	6µl	R57028	N/A